

Operating Manual

Fluent®



Title:	Fluent Operating Manual			Part number:	30134933.13
ID:	399706, en, V2.2			Translated from:	n.a.
Version:	Revision:	Issue:	Document History:		
1	0	2017-08-30	First Edition		
1	1	2017-09-18	2.3 Application Risks		
1	2	2017-11-07	1.3 Intended Use		
1	3	2018-03-05	1.3 Intended Use		
1	4	2018-04-16	3.3 Dimensions and Weights		
1	5	2018-07-20	4.6.2 FCA Gripper		
1	6	2019-02-26	2.2 General Safety Information		
1	7	2020-10-22	Information on Mix & Pierce added in sections 1.8.3, 2.5, 2.8.1, 4.6, 4.6.3, 6.6, 7.4.1, 7.4.2, 7.5.4, 8.2.3, 8.3.13–15 Information on Tube Rotator added in sections 4.6.3, 6.5.3, 7.4.1, 7.4.2, 7.5.5, 7.5.6, 8.2.3, 8.3.14 Sections 4.6.4 and 7.5.24 added. Sections 2.5, 2.9, 2.10, 4.2.2, 4.6, 6.3.5, 7.1, 7.4, 7.5.16, 8.2.4, and 8.3.13–20 updated		
1	8	2020-12-01	Sections 2.2 and 2.8 updated		
1	9	2021-02-15	Sections 2.7.2, 7.1, 7.2, 7.4.2, 7.5.15, 7.5.24, and 8.2.3 updated Sections 2.7.5, 4.6.4, and 7.5.25 added		
1	10	2021-11-25	Section 3.7 updated Section 6.5.1 updated		
2	0	2022-03-02	Section 3.1 updated		
2	1	2022-04-14	Sections 5.2.5 and 6.6 added Section 6.5.1 updated		
2	2	2022-11-16	Section 4.6.5 added Section 2.2 and 3.4 updated		

© 2022, Tecan Trading AG, Switzerland, all rights reserved.

Information contained in this document is subject to change without notice.

Table of Contents

1 About This Manual	11
1.1 Scope of This Manual	11
1.2 Manufacturer.....	11
1.3 Intended Use	11
1.4 Area of Use/ Area of Application.....	12
1.5 Improper Use	12
1.6 Warranty	12
1.7 Trademarks.....	12
1.8 Reference Documents	12
1.8.1 Instrument Manuals	13
1.8.2 Software Manuals	13
1.8.3 QC Kit Manuals.....	13
1.8.4 Other Reference Documents	13
1.9 Compliance with Laws and Standards.....	13
1.10 Document Conventions	13
2 Safety.....	14
2.1 Safety Message Conventions	14
2.1.1 Signal Words.....	14
2.1.2 Safety Symbols	15
2.2 General Safety Information.....	16
2.3 Camera Privacy Statement.....	17
2.4 Application Risks	18
2.5 Operating Company.....	22
2.6 Method and Process Validation	23
2.7 User Qualification	24
2.7.1 Operator	24
2.7.2 Key Operator.....	24
2.8 Safety Elements.....	25
2.8.1 Safety Panels.....	25

2.8.2	Safety Panel Sensors/Active Stop	30
2.8.3	Instrument Door Locks (optional)	30
2.8.4	Cabinet Door Locks	30
2.8.5	External Door Locks	30
2.9	Product Safety Signs	31
2.9.1	Mix & Pierce Workstation	35
2.10	Laser Radiation	36
2.10.1	Laser Radiation Devices	36
2.11	Optical Radiation (UVC)	37
2.12	Decontamination Declaration	38
3	Technical Data	39
3.1	Type Plate	39
3.2	Serial Number Label	40
3.3	Dimensions and Weights	41
3.4	Power Supply	42
3.5	Data and Power Connections	43
3.6	Environmental Conditions	44
3.7	Emission and Immunity	45
4	Description of Function	46
4.1	Overview	46
4.2	Deck	46
4.2.1	Carriers	47
4.2.2	Deck Trays	47
4.2.3	Placement on Segment	48
4.2.4	Segment Position	49
4.3	Robotic Arms	49
4.3.1	Flexible Channel Arm (FCA)	50
4.3.2	Multiple Channel Arm (MCA)	51
4.3.3	Robotic Gripper Arm (RGA)	51
4.3.4	Arm Accessories	52

4.4	Liquid System (Liquid FCA)	53
4.5	Wash System (MCA)	54
4.6	Options and Devices	54
4.6.1	Fluent ID Tube Barcode Scanner	55
4.6.2	FCA Gripper	56
4.6.3	Mix & Pierce	58
4.6.4	Frida Reader	61
4.6.5	Phase Separator	62
5	Control Elements	63
5.1	Operating Elements	63
5.2	User Interface	64
5.2.1	Navigation Path	64
5.2.2	Working Area	64
5.2.3	Display, Option and Action Buttons	65
5.2.4	Method Recovery Buttons	66
5.2.5	DeckCheck Buttons	67
5.3	Error Signals and Instrument Status	69
5.4	Fluent ID Status LEDs	71
6	Operation	72
6.1	Safety Instructions for This Chapter	72
6.2	Operating Modes	72
6.3	Putting into Operation	73
6.3.1	Switching On the Instrument	73
6.3.2	Starting FluentControl	74
6.3.3	User Login	74
6.3.4	Placing Segments	75
6.3.5	Removing Segments	76
6.3.6	Loading Standard Runners	78
6.3.7	Checking the Deck Layout	80
6.4	Before Starting a Method	81

6.4.1	Checking the Tubing on System Liquid Container and Waste Container	84
6.4.2	Checking Waste Container Tubing	85
6.5	Running a Method	85
6.5.1	Starting a Method.....	86
6.5.2	Loading and Unloading Fluent ID Runners.....	89
6.5.3	Loading and Unloading of Tube Rotator Runners	93
6.5.4	Resetting Errors	96
6.6	DeckCheck Operation.....	97
6.7	Method Recovery.....	99
6.7.1	Switching to Method Recovery Mode	100
6.7.2	Recovering a Method Run	100
6.8	Switching Off the Instrument.....	101
7	System Care.....	103
7.1	Decontamination.....	103
7.2	Cleaning Agents	104
7.2.1	Cleaning Agents Specifications	104
7.2.2	Commercial Cleaning Agents	104
7.3	System Care Mode	105
7.3.1	Switching to System Care Mode.....	105
7.3.2	Resetting Errors	106
7.4	System Care Tables	106
7.4.1	Daily System Care	107
7.4.2	Weekly System Care	111
7.4.3	Monthly System Care.....	113
7.4.4	Periodic System Care	113
7.4.5	Yearly System Care	114
7.4.6	Biennial System Care	114
7.5	System Care Activities.....	114
7.5.1	Moving the Instrument on a Cabinet within the Laboratory.....	114
7.5.2	Cleaning Disposable Tip Cone	117
7.5.3	Cleaning Fixed Tips	118

7.5.4	Cleaning Piercing Tips	118
7.5.5	Cleaning the Tube Rotator	118
7.5.6	Cleaning the Tube Rotator Wash Station	121
7.5.7	Cleaning Runners and Segments	123
7.5.8	Cleaning Deck Trays.....	123
7.5.9	Replacing Fluent ID Reflector Foil	125
7.5.10	Applying Fluent ID Reflector Foil on DiTi Waste Slide.....	126
7.5.11	Cleaning Safety Panels.....	127
7.5.12	Cleaning Disposable Tip Waste and Wash Station Unit	127
7.5.13	Cleaning Disposable Tip Waste Slide.....	128
7.5.14	Changing Disposable Tip Waste Bag	130
7.5.15	Cleaning Liquid Path.....	131
7.5.16	Connecting the System Liquid Container and Waste Container.....	132
7.5.17	Cleaning the System Liquid Container and Waste Container.....	133
7.5.18	Checking Tightness of Syringes	133
7.5.19	Checking Gaskets (MCA)	134
7.5.20	Replacing Gaskets (MCA)	135
7.5.21	Cleaning Fixed Tip Block (MCA).....	138
7.5.22	Cleaning Plate Adapter (MCA).....	138
7.5.23	Cleaning Arm Guide.....	138
7.5.24	Tightening a DiTi Cone	138
7.5.25	Frida Reader	139
8	Troubleshooting	141
8.1	Safety Instructions for This Chapter	141
8.2	Troubleshooting Tables	141
8.2.1	Instrument Troubleshooting	141
8.2.2	Flexible Channel Arm (FCA) Troubleshooting	142
8.2.3	Mix and Pierce	143
8.2.4	Multiple Channel Arm (MCA) Troubleshooting	146
8.2.5	Robotic Gripper Arm (RGA) Troubleshooting	148
8.2.6	Wash System Troubleshooting.....	149

8.2.7	Fluent ID Troubleshooting.....	150
8.2.8	Software Troubleshooting	151
8.3	Troubleshooting Activities.....	152
8.3.1	Connecting Wash Station (MCA).....	152
8.3.2	Position Labware	153
8.3.3	Checking Inline Filter (Air FCA)	154
8.3.4	Changing Inline Filter (Air FCA).....	154
8.3.5	Removing DiTi Cone (Air FCA).....	156
8.3.6	Assembling DiTi Ejector Tube (Air FCA).....	157
8.3.7	Installing DiTi Cone (Air FCA).....	158
8.3.8	Removing the DiTi Option (FCA)	160
8.3.9	Installing the DiTi Option (FCA)	160
8.3.10	Checking Fixed Tips	161
8.3.11	Removing Fixed Tips	162
8.3.12	Installing Fixed Tips	163
8.3.13	Removing Piercing Tips	165
8.3.14	Installing Piercing Tips	168
8.3.15	Retracting Stuck Piercing Tips.....	172
8.3.16	Checking Gripper Finger Alignment.....	173
8.3.17	Basic Gripper Fingers Alignment for FES Gripper Fingers	174
8.3.18	Basic Gripper Fingers Alignment for Fixed Gripper Fingers	175
8.3.19	Advanced Gripper Fingers Alignment for FES Gripper Fingers.....	176
8.3.20	Advanced Gripper Fingers Alignment for Fixed Gripper Fingers.....	176
8.3.21	Checking Segment.....	178
8.3.22	Removing Positioning Pins	179
8.3.23	Replacing Lock Pins and Positioning Pins.....	180
9	Packing, Unpacking, Transport, Storage and Disposal.....	182
9.1	Packaging Labels	182
9.2	Disposal	183
9.2.1	Local Requirements European Union	183
9.2.2	Local Requirements People's Republic of China.....	183

9.2.3 Other Requirements.....	184
10 Customer Support.....	185
10.1 Contacts.....	185
Abbreviations.....	188

1 About This Manual

This Operating Manual describes the Fluent and provides all the information required for operating it safely and for maintaining it in good working order. This manual must be read carefully before performing any work on the Fluent and before using it.

This chapter outlines the purpose of this manual and specifies the product referred to. Furthermore, it explains the use of symbols and conventions as well as further general information.



This Operating Manual contains no software description. For more information on the software please consult the corresponding software manual. Refer to section [“Reference Documents”](#) [► 12].

1.1 Scope of This Manual

This manual applies to:

- Fluent 480 (part number 30042011)
- Fluent 780 (part number 30042021)
- Fluent 1080 (part number 30042031)

1.2 Manufacturer

Address of
Manufacturer



Tecan Schweiz AG
Seestrasse 103
CH-8708 Männedorf
Switzerland

1.3 Intended Use

The Fluent is a fully automated laboratory liquid-handling platform for general laboratory use. It is intended for routine laboratory tasks, such as pipetting, liquid handling, and robotic manipulation of labware in certain test procedures.



The Fluent Gx Assurance Software is required for this intended use.

Certain options from Tecan and third-party devices that can be used with Fluent are for research use only (RUO).

If a research-use-only option or device is integrated with Fluent or if the Fluent Gx Assurance Software is not installed, the intended use changes to:

The Fluent is a fully automated laboratory liquid-handling platform for research and industrial applications. It is intended for routine laboratory tasks, such as pipetting, liquid handling, and robotic manipulation of labware in certain test procedures (e.g., cell-based assays, biochemical assays and compound management). This instrument is not intended for use in clinical diagnostics procedures.

1.4 Area of Use/ Area of Application

Fluent may be used in a variety of laboratory environments according to its intended use.

In each environment the individual laboratory is responsible for the validation of the Fluent instrument together with the specific liquids and labware used in the laboratory's application workflow or method.

1.5 Improper Use

Improper use may prejudice the Fluent safety concept.

- The Fluent must not be used with options or components that are not approved by Tecan.
- The Fluent is not explosion-proof and should not be installed in locations where there is a hazard of explosion.
- The Fluent should not be used in the absence of functional safety devices.

1.6 Warranty

The Fluent must not be used with components that are not approved by Tecan.

The use of unapproved components may impair the safety concept of the Fluent.

The use of unapproved components would invalidate any warranty of safety and compliance to national and international standards, as required for UL/CSA certification, by EC directives, etc.

1.7 Trademarks

The product names, whether registered or unregistered trademarks, mentioned in this manual are reproduced solely for identification purposes and remain the exclusive property of their respective owners. For simplicity reasons, the trademark symbols such as ® and ™ are not repeated in the manual.

1.8 Reference Documents

This section provides a list of the documents that are needed or may be useful when using the Fluent.

The Doc IDs listed below are root numbers. Therefore, they do not contain information about the language, document version, or the medium (data storage medium, hard copy, downloadable file, etc.) of the document.



On the basis of your order configuration, the Operating Manuals for optional equipment apply as well.

Check the scope of the corresponding document to ensure that you are in possession of the correct version.

The Doc ID does not refer to ordering information. When placing orders, please refer to the number on the binder, CD casing, etc.

1.8.1 Instrument Manuals

- Fluent® Operating Manual (Doc ID 399706)
- Fluent® Reference Manual (Doc ID 399937)

1.8.2 Software Manuals

- Tecan Sample Tracking Add-on Software Manual (Doc ID 393933)

1.8.3 QC Kit Manuals

- QC Kit Application Manual (Doc ID 397069)
- QC Kit Application Software Manual (Doc ID 397070)

1.8.4 Other Reference Documents

- Fluent® Carousel Operating Manual (Doc ID 398350)
- HEPA hood (Doc ID Bigneat 70072)
- Frida Reader™ Application Manual (Doc ID 401882)
- Te-Shake™ Operating Manual (Doc ID 391496)
- Te-VacS™ Operating Manual (Doc ID 391236)
- Fluent® Stacker Operating Manual (Doc ID 398658)
- MIO2 Operating Manual (Doc ID 394934)

1.9 Compliance with Laws and Standards

The following declarations and certifications apply to Fluent:

- EC Declaration of Conformity with applicable EU Directives (CE mark)
- Canadian Standard Association Certification (CSA mark)
- (IECEE) CB Scheme Certification (CB mark)

For more detailed information about the marking, refer to section Type Plate.

1.10 Document Conventions

Cross-References

Cross-references appear as follows—e.g.:

Refer to section “[Safety](#)” [▶ 14]

- “Safety” refers to the corresponding section header
- The page number is given in square brackets

Prerequisites

Prerequisites appear as follows—e.g.:

- ✓ “General Safety Information” has been read.

Tips

Additional tips appear as follows—e.g.:



For safety conventions and symbols refer to chapter “[Safety](#)” [▶ 14].

Illustrations

The illustrations may show component versions which are not relevant to your Fluent.

2 Safety

This chapter describes the safety concept of Fluent, provides general rules of correct behavior, and warnings concerning hazards associated with the use of the Fluent.

2.1 Safety Message Conventions

2.1.1 Signal Words

Tab. 1: Signal Words

Signal Word	Meaning
 DANGER	Indicates a hazardous situation that, if not avoided, will result in death or serious injury.
 WARNING	Indicates a hazardous situation that, if not avoided, could result in death or serious injury.
 CAUTION	Indicates a hazardous situation that, if not avoided, could result in minor or moderate injury.
 NOTICE	Indicates a situation that is not hazard-related but, if not avoided, could result in damage to or malfunctioning of the equipment, or incorrect process results.

2.1.2 Safety Symbols



Crushing of hands



General warning



Laser beam



Optical radiation



Biohazard



No heavy load

2.2 General Safety Information

WARNING

Fluent is designed and built in accordance with the present state-of-the-art technology and the recognized technical safety regulations. Nevertheless, risks to users, property and the environment can arise if the Fluent is used without due care and attention.

The safety of all users and personnel depends on the strict observation of these safety instructions and awareness of the safety-related warnings provided in this manual.

- Please pay great attention to the following general safety information.
 - This manual must always be available to all persons performing the tasks described herein.
-
- Legal regulations, such as local, state and federal laws concerning the use or application, as well as the handling, of dangerous materials in connection with the Fluent must be strictly followed.
 - The operating company is responsible for defining instructions in accordance with company procedures and local legal requirements. The instructions provided by the operating company must be strictly observed.
 - Observe the correct environmental conditions for storage and operation.
 - Structural changes to the safety devices are forbidden.
 - Damaged safety devices must be replaced immediately as described in this manual.
 - The Fluent must not be modified in any way without prior consultation and written approval of Tecan. Authorized modifications to the system may only be performed by an FSE certified for the repair and upgrading of the Fluent. Tecan will reject any claim resulting from unauthorized modifications.
 - Fire hazard caused by the improper use of the Fluent. The Fluent should not be installed in locations where there is a hazard of explosion.
 - Chemical, biological, and radioactive hazards can be associated with the substances used or the samples and reagents processed with the Fluent (e.g., during loading and unloading). The same applies to waste disposal.
 - Always be aware of possible hazards associated with these substances.
 - Use appropriate protective clothing, safety goggles, respirators, and gloves.
 - The handling of substances and the disposal of waste may be subject to local, state, or federal law, or to regulations with regard to health, environment, or safety. Strictly observe the corresponding provisions.
 - Any contamination must be dealt with immediately as described in this manual.
 - The user is responsible for ensuring that the Fluent is always operated under proper conditions, and that maintenance, service, and repair tasks are performed with care, on schedule, and only by authorized personnel.
 - Risk of incorrect measuring results. After system care or maintenance has been performed, operation must only be resumed after the correct system operating conditions have been verified.

- Always use recommended consumables within expiration date and original spare parts for maintenance and repair to assure good system performance and reliability.
- Injury could result if skin comes in contact with the instrument or system liquid.
 - Always wear protective clothing according to GLP.
- Heavy load! Do not lift the instrument.
- Fire hazard caused by flammable liquids or system liquid.
 - Avoid the formation and accumulation of flammable vapors.
- Do not operate the system without deck trays and deck segments.
- Deck trays capture liquid spills that may occur in the manual deck loading area. The system should be operated with as many deck trays as possible installed below the deck to collect all liquid spills. Do not operate the system without deck trays.
- If carry-over is not tolerated, the use of disposable tips with filters is strongly recommended.
- Possible Crash. Do not place devices without Tecan model data on the deck.
- Extension 300 is designed for a maximum load of 40 kg (88 lbs.), and only for use with options that are lighter than 40 kg (88 lbs.).
- The Fluent is supplied with a biohazard safety sign which should be applied by the user in the event of use of biohazardous substances. Apply the label on the front door in a position visible to the user and convenient for the application. Refer to section Product Safety Signs.
- The ethernet cable of the DeckCheck cameras will be installed by an FSE and must be installed on the Fluent PC at all times (EMC). The Ethernet interface is not allowed to be connected to a network.
- For California residents only: This product can expose you to chemicals such as lead which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov/product.

2.3 Camera Privacy Statement

The Fluent system is equipped with cameras mounted on the inside front profile. The cameras are focused on the deck and rear deck. Views downwards through the acrylic-glass side panels are possible.

- The user is responsible for advising people in the room that cameras are in operation.
- The user is responsible for ensuring that the personnel could not be identified from pictures taken, for example if the instrument is adjacent (side-on) to a desk space or if rear or side panel cut-outs are made or if an acrylic-glass panel is used to replace the rear wall.

2.4 Application Risks

System function / Module	Possible failure mode	Potential failure effect	Possible / Potential cause	Labeling or mitigation
System	Insufficient maintenance	Safety or health of users: Potential contamination of instrument	Use error: Neglecting Operating Manual or maintenance instructions	The Operating Manual informs the user about the use of appropriate consumables and preventive maintenance instructions The Operating Manual contains a reference that user shall wear protective clothing, gloves and goggles according to GLP
System	Fire	Safety or health of users: Fire in operators lab (instrument burning)	Gas from volatile flammable liquids; spark from electronic board spreading	The Operating Manual includes the following information: The instrument is not explosion safe and the customer shall ensure that there is no high vapor concentration.
Module FCA & Air FCA	Wear out of Z-axis mechanics (Above average usage)	Safety or clinical conditions of sample: Potentially wrong Z-positioning in labware	Above average usage of device in combination with usage of disposable tips High percentage of piercing steps in application	The system informs the user if they have reached 90% of the expected lifetime of the axis.
Module FCA & Air FCA	Abrasion of X-drive cog-wheel (Above average usage)	Safety or clinical conditions of sample: Potential contamination of samples with polyamide particles	Above average usage of device in combination with placing of critical labware on rear of instrument	Avoid placing particle-sensitive elements (e.g., samples and reagents) on rear of the instrument or place particle protection on top of labware (i.e., lids)
Module FCA & Air FCA	Interfering signals due to septum piercing	Safety or clinical conditions of patient sample: wrong cLLD leading to air aspiration and potentially false results	Interaction of the tip with the septum / foil	The Reference Manual informs the user to work only with non-conductive foils for piercing applications in conjunction with liquid level detection on the FCA and Air FCA The Reference Manual informs the user to validate liquid detection in combination with piercing for FCA and Air FCA

System function / Module	Possible failure mode	Potential failure effect	Possible / Potential cause	Labeling or mitigation
Module FCA & Air FCA	Wrong sample treatment, wrong cLLD due to foam or bubbles in reagent vial	Safety of process: wrongly processed samples	Bubbles or foam in the reagent vial cause a wrong cLLD and potential aspiration of air with FCA or Air FCA	The Reference Manual informs the user to validate the application / process.
Module FCA & Air FCA	Tip blockade	Safety or clinical conditions of patient sample: Potentially wrong pipetted volume	Aspiration at bottom of well (blockade of tip)	The Operating Manual informs the user to validate the application to prevent aspiration too close to Z-max level of custom labware
Module FCA & Air FCA	FCA tubing system: Growth of microorganisms	Safety or clinical condition of patient sample: wrong pipetted volume or contamination of samples	Growth of microorganisms (biofilm on the inner surface)	The Operating Manual informs the user to use deionized water as system liquid for FCA and recommends daily maintenance to flush the system (also list of allowed cleaning reagents)
Module MCA	Overflow of sample liquid in microplate during pipetting	Safety or clinical condition of patient sample: Potential cross-contamination of samples (overflow)	Wrongly defined Z-levels by user (e.g., aspiration from Z-max position)	The Reference Manual contains instructions for the definition of safe positions for aspiration and dispensation
Module MCA	Samples completely or partially miss the intended position in free dispense mode	Safety or clinical condition of patient sample: Potential cross-contamination	Electro-static charges at the tip end due to instrument usage outside specified conditions lead to sample remaining hanging on the tip or to uncontrolled spraying	The software only contains default Liquid Classes with correctly defined Z-dispense levels (inside well) The Operating Manual contains specified operation conditions for MCA Liquid Handling especially in instruction about min. required humidity The Reference Manual contains a warning to set the dispense height to be inside of the well

System function / Module	Possible failure mode	Potential failure effect	Possible / Potential cause	Labeling or mitigation
Module MCA	Mixing air instead of liquid (sample / reagent) for Mix-pipetting	Safety or clinical condition of patient sample: Samples potentially incorrectly processed resulting in false results	Inappropriate tracking parameters due to wrong combination of tips and microplates	The Software helps the user to compare the real and virtual worktable by showing the name of the labware in the virtual worktable. The Mechanical design defines a unique color design (Tip type specific) and Labeling (for Filter and Non Filter) of DiTi Boxes The Operating Manual contains instructions to check the worktable layout before starting a process
Disposable Tip specific	Get DiTis: Incorrect tip type mounted	Safety or clinical condition of patient sample: Potentially no or short sample aspirated Potential cross-contamination of samples	Use error: Incorrect deck layout: user puts tip box at wrong position: tips are shorter than expected Incorrect deck layout: user puts tip box with unfiltered tips instead of filtered tips on worktable Incorrect deck layout: user puts tip box at wrong position: tip has smaller volume than expected (e.g., 100 µl instead of 200 µl); tip length as expected; liquid aspirated in MCH	The software helps the user to compare the real and virtual worktable by showing the name of the labware in the virtual worktable. The mechanical design defines a unique color design (tip type specific) and labeling (for filter and non filter) of DiTi boxes The Operating Manual contains instructions to check the worktable layout before starting a process The mechanical design ensures visibility of the white filter The Reference Manual contains information about color coding of DiTi boxes and difference of length, and filtered DiTis
Disposable Tip specific	Incomplete dropping of tips: Some contaminated tips remain hanging to the head and fall onto sample plates	Safety or clinical condition of patient sample: Potential cross-contamination	Electro-static charges caused	The Operating Manual contains specified operation conditions for MCA liquid handling—especially instruction about min. required humidity The Operating Manual contains a reference that disposable tips are not intended for reuse

System function / Module	Possible failure mode	Potential failure effect	Possible / Potential cause	Labeling or mitigation
Module RGA	Plate loss due to crash with misaligned labware	Safety of process: Plate loss, loss of samples	If more than 4 microplates are stacked, misalignment can occur during transport	The Application Software Manual informs the user that plate movements have to be validated.
Module FluentControl Software	Worktable-Base: wrong DiTi status reported	Safety of Process: Cross-contamination / Wrong results	Cross-contamination due to wrong information about usage status of tips	Do not use "Set Tips Back" if failure mode leads to high severity risk
Module FluentControl Software	Core.Scripting.Programming Set-Variable at run time: wrong value	Process safety: wrong results	Error in software: variable is set to wrong value	Validate the application for the specific variable source, destination and ranges
Module FluentControl Software	Core.Scripting.Programming Query-Variable at run time or script start: wrong UI presentation / acceptance of UI value	Process safety: wrong results	Numeric value is formatted or converted wrongly in UI	Validate the application for the specific variable source, destination and ranges
Module FluentControl Software	Core.Scripting.Programming Import-Variable at runtime: wrong value imported	Process safety: wrong results	Wrong value is retrieved from import source	Validate the application for the specific variable source, destination and ranges
Module FluentControl Software	Core.Scripting.Programming Export-Variable at run time: wrong value exported to file	Process safety: wrong results	Wrong value is written to export file	Validate the application for the specific variable source, destination and ranges

System function / Module	Possible failure mode	Potential failure effect	Possible / Potential cause	Labeling or mitigation
Module FluentControl Software	API: Get/set variable or resolve expression fails	Process safety: wrong results	Wrong variable value retrieved or wrong value assigned / wrong expression result returned	Validate the application for the specific variable source, destination and ranges
UVC light	Incorrect usage in application	Lack of effectiveness	Incorrect usage in application	Refer to specific instructions in section “Optical Radiation (UVC)” [▶ 37]
Tube rotator/piercing tips (Mix & Pierce)	Incorrect usage in application	Lack of effectiveness	Incorrect usage in application	Refer to specific instructions in section “Mix & Pierce” [▶ 58]
Frida Reader	Incorrect usage in application	Lack of effectiveness	Incorrect usage in application	Refer to specific instructions in section “Frida Reader” [▶ 61]
Any	Ineffective usage in application	Lack of effectiveness in application	Lack of system Care	Refer to part specific instructions in chapter System Care
Processing potentially hazardous materials	Contamination with potentially hazardous materials	Potential risks to users, property and the environment	Lack of adherence to General safety information	Refer to part specific instructions in section General Safety Information

2.5 Operating Company

The operating company must ensure that the Fluent and in particular the safety features, function properly and that all the personnel in contact with the instrument are adequately trained.

Responsibilities

- Method and process validation.
- Defining the processes in compliance with the Standard Operating Procedures.
- Ensuring that installation and operational qualifications (IQ OQs) have been completed.
- Ensuring that all personnel in contact with the Fluent are adequately trained.
- Ensuring the availability of appropriate protective clothing and equipment.
- Ensuring the maintenance and safe operation of the Fluent.
- Requiring adherence to laboratory safety regulations and directives.

2.6 Method and Process Validation

While performing method and process validation, pay attention to the following:

- If using fixed tips with MCA or FCA, ensure that the wash procedure is effective for the expected sample concentration range and assay sensitivity.
- Check that pipetted volumes meet the precision and accuracy requirements of the process being automated.
- When using non-Tecan or custom labware and aspirating with tracking, ensure that the container definition is correct (i.e., the appropriate speed is used for tracking) to avoid air aspiration.
- The Phase Separator functionality has been verified for use with standard Tecan 1 ml disposable tips and 1 ml Tecan wide-bore disposable tips. For more information on supported Tecan consumables, refer to the Reference Manual (see [“Reference Documents” \[▶ 12\]](#)).
- Validate liquid detection on the Fluent Stacker transfer station.
- Validate the correct usage of the MCA wash station by the application.
- Validate the application with regard to correct pipetting volumes and tracking.
- Validate the application to prevent aspiration too close to Z-max of custom labware.
- Validate piercing applications with regard to the downholders needed (active or passive).
- If chemicals and labware are not removed, the impact of UVC light on chemicals and labware present on the deck has to be evaluated and the assay validated.
- Include a manual post-run check for correct pipetting volumes.
- Personnel must be informed regarding the camera privacy statement (refer to [“Camera Privacy Statement” \[▶ 17\]](#)).

**Responsibility
of the Key
Operator**

2.7 User Qualification

The laboratory personnel must be fully qualified and trained to operate the Fluent. The work described in this Operating Manual must only be performed by authorized personnel with the qualifications prescribed below.

Laboratory personnel must:

- have suitable technical training,
- be familiar with the laboratory safety regulations and directives,
- be familiar with the instructions for the safety elements of the instrument,
- use protective clothing and equipment,
- be familiar with and adhere to good laboratory practices,
- and have read and understood the instructions in the Operating Manual.

Tecan recommends that the operator attends an operator training course. Please ask the Tecan Customer Service about available courses. Refer to section [“Customer Support”](#) [► 185].

2.7.1 Operator

The operator (lab technician) works for the operating company.

Required Skills

- No specific application or system knowledge
- Command of local languages
- Command of English is preferable

The operator has application software access rights allowing him to run methods and perform system care.

2.7.2 Key Operator

The key operator (application specialist) supports the operating company or works for the same company.

Required Skills

- Extensive application knowledge
- Limited system knowledge
- Command of local languages
- Command of English
- In-depth knowledge of the corresponding software manual

Responsibilities

- Instructing the operator
- Writing, running and validating methods
- Helping the operator to solve problems with the instrument

2.8 Safety Elements

CAUTION

Moving parts

The protection and safety elements installed on the Fluent must not be removed, disabled or overridden during operation.

- If any devices are removed (e.g., for maintenance work), all protection and safety devices must be reinstalled, re-enabled and checked before resuming operations.

Safety panels and safety sensors are integral parts of the Fluent, whereas instrument door locks and cabinet door locks can be included only in certain system configurations.

2.8.1 Safety Panels

Fluent is protected with safety panels:

The **front safety panel** can be opened and is fitted with door sensors that trigger an active stop. The front safety panel can be locked with optional door locks.

The **diluter panel** can be opened without affecting Fluent operation.

The **top and side safety panels** are fixed.

2.8.1.1 Front Safety Panels

The front safety panel prevents direct access to the robotic arms and to the elements on the instrument deck during operation. This is for the benefit of personal safety and improves method security. In addition, the front safety panel protects the user against spilling sample or reagent. There are different types of front safety panels.

Full Front Safety Panel




Fig. 1: Full Front Safety Panel

The full front safety panel has the following features:

- No access to moving parts (moving parts, mechanical hazards)
- Protection of the samples against outside influence (method safety)
- Protection against spilling sample or reagent



With full front safety panels, only batch-wise loading is possible.

**Full Front
Safety Panel
(UVC)**




Fig. 2: Full Front Safety Panel (UVC)

The full front safety panel (UVC) has the following features:

- No access to moving parts (moving parts, mechanical hazards)
- Protection of the samples against outside influence (method safety)
- Protection against spilling sample or reagent
- Protection against optical radiation (UVC)



With full front safety panels, only batch-wise loading is possible.

⚠ CAUTION

Moving parts!

Moving MCA, FCA and Air FCA can cause hand injuries when reaching through the half front safety panel or of the front safety panel with expansion into the instrument during a run.

- Do not reach into the instrument during a run.
-

Half Front Safety Panel




Fig. 3: Half Front Safety Panel

The half front safety panel has the following features:

- Restricted access to moving parts (moving parts, mechanical hazards)
- Protection against spilling sample or reagent



With the half front safety panel the operator has restricted access to the instrument deck. Loading and unloading runners is possible without opening the panel; i.e., the operator is enabled to reload samples or reagents during the method run.

Front Safety Panel with Expansion




Fig. 4: Front Safety Panel with Expansion

The front safety panel has the following features:

- Restricted access to moving parts (moving parts, mechanical hazards)
- Protection against spilling sample or reagent
- Allows for the use of a front DiTi waste station, which protrudes from the deck and requires a downward-facing opening in the front safety panel.



With the front safety panel with expansion, only batch-wise loading is possible.

2.8.1.2 Safety Panels for Optional Devices

If an optional device is added to, or removed from, the side of the Fluent, an appropriate side safety panel must be installed. Please consult the [“Customer Support” \[185\]](#).

2.8.2 Safety Panel Sensors/Active Stop

The Fluent safety concept assumes that the front safety panel is always closed when the instrument is running.

Active Stop

As soon as the front safety panel is opened, an active stop is triggered by door sensors. This means that all arm movements come to a halt for safety reasons. The operator must close the safety panel and resume the program to continue the process. **CAUTION! Unauthorized modifications to door sensors are prohibited.**

The entire run completes “with warnings”. It is recommended to check the errors and warnings before releasing the run if the run does not complete successfully.



The following devices will not be interrupted by an active stop: Tecan Incubator, Magellan, Te-Shake, Fluent Stacker. Interruption of other devices will depend on the device driver.

2.8.3 Instrument Door Locks (optional)

Two optional door locks can prevent the front safety panel from being opened and protect the ongoing process. This prevents unwarranted interruption of the process run. To stop a process, a pause request can be entered by means of the touchscreen.

2.8.4 Cabinet Door Locks

If an RGA long axis has access below the deck, the cabinet door closest to the access point must be equipped with a door lock sensor option. If more than one access point below the deck is implemented or if the access point is changed during the life of the instrument, then each door near the access point must be equipped with a door lock sensor.

If the instrument has a HEPA hood, all cabinet doors must be equipped with a door lock sensor.

2.8.5 External Door Locks

External door locks will be implemented on Fluent installations in an external enclosure. The door panels of the external enclosure replace the mechanical safety function of the Fluent front safety panel and the cabinet doors, and the external door docks with integral sensors replace the door sensor and door lock functions of the Fluent front safety panel and cabinet doors.



External door locks do not allow an ActiveStop. To stop or pause the process a pause request can be entered by means of the touch screen.

2.9 Product Safety Signs

Safety signs are affixed to the Fluent for safety purposes. Damaged, lost or illegible safety signs must be replaced immediately as illustrated. For the meaning of safety symbols refer to section “[Safety Message Conventions](#)” [14].

Standard Instrument




Fig. 5: Standard instrument

UVC




Fig. 6: Instrument with UVC




Fig. 7: Inside view

Biohazard




Fig. 8: Biohazard



The Fluent is supplied with a biohazard safety sign which should be applied by the user in the event of use of biohazardous substances.

Apply the label on the front door in a position visible to the user and convenient for the application.

**Instrument with
Half Front
Safety Panel**




Fig. 9: Instrument with Half Front Safety Panel

**Instrument with
Front Safety
Panel with
Expansion**




Fig. 10: Instrument with Front Safety Panel with Expansion

MCA



Fig. 11: Safety sign on MCA

Deck Extension




Fig. 12: Deck extension

2.9.1 Mix & Pierce Workstation

FCA Safety Shield



Fig. 13: Safety shield

2.10 Laser Radiation

Fluent can be equipped with laser barcode scanners. The laser radiation from these barcode scanners is a low-power, collimated beam in the visible spectrum. The laser classes of each barcode scanner, and of the entire Fluent system, are indicated on the laser safety label affixed to the corresponding hardware.

All modules with lasers are marked with the appropriate laser safety labels.

The Fluent instrument has been tested and certified according to IEC 60825-1:2007 and IEC 60825-1:2014.



CAUTION

The Fluent is a class 1 laser product pursuant to IEC 60825-1:2014 that emits laser radiation.

Dazzle, flash-blindness and afterimages may be caused by the laser beam.



- Do not stare into the laser beam or into its reflections.


2.10.1 Laser Radiation Devices

A stand-alone barcode scanner can be mounted on a device.

Please ensure that the safety label is correctly affixed to the barcode scanner at all times:

- Explanatory Laser Radiation Label (A): Identifies a CLASS 2 LASER PRODUCT according to IEC 60825-1 that contains an embedded visible low power laser barcode scanner. Instructs the user not to stare into laser beam or its reflection.
- Class 2 lasers are only operated when the system is running and have no interface to the operator.

Label placement	Explanation
 	Stand-alone barcode scanner mounted on a hotel: label located underneath the scanner.
 	Stand-alone barcode scanner mounted on Robotic Gripper Arm: label located on the scanner.

Label placement	Explanation
	Fluent ID: label located on the rear side of the scanner housing.
	Fluent ID: label located on the side of the scanner housing.

2.11 Optical Radiation (UVC)

The Fluent can be equipped with an optional HEPA hood which includes a UVC light, or a separate UVC light option.

Exposure to UVC light radiation must be avoided as it can lead to injury. The UVC light switches off automatically when the front safety panel is opened, and in the case of the UVC light option also when the diluter cover is opened. Special UVC-resistant safety panels are installed on the Fluent in conjunction with UVC light.

UVC light can be used in decontamination procedures. The suitability and effectiveness of using UVC for individual processes must be validated by the user.



Please also refer to the manual provided by the HEPA hood manufacturer.

2.12 Decontamination Declaration

In addition to regular system care, and in accordance with standard laboratory regulations, the Fluent and its parts and accessories must be thoroughly decontaminated in the following circumstances:

- Before any maintenance or service work is performed on the Fluent and, in particular, before an FSE intervention on the Fluent
- In the event of accidents (e.g., crash, spillage, etc.)
- Before returning the Fluent or its parts or accessories, to Tecan (e.g., for repair)
- Prior to storage
- Prior to disposal
- In general, before moving the Fluent or its parts from its location

The owner of the instrument has full responsibility for the effective decontamination of all the equipment.

Before any intervention on the Fluent by an FSE, and before returning the Fluent or its parts or accessories to Tecan, the owner of the instrument must complete and sign the Decontamination Declaration form, confirming that the decontamination has been performed in accordance with good laboratory practice guidelines. Contact your local service organization to obtain this form and refer to section Decontamination.



Tecan reserves the right to refuse to deal with any Fluent or its parts or accessories that is not accompanied by the Decontamination Declaration form.

3 Technical Data

3.1 Type Plate




Fig. 14: Type plate

The type plate is on the rear side of the Fluent and contains the following information:

Identification data	Model
	REF: Ordering information (material number and revision level)
	Date of manufacture (YYYYMMDD)
	SN: Serial number
Technical data	U, f: Supply voltage (Volts), frequency (Hertz)
	P: Power consumption (W)
Address data	Manufacturer's name and address
Conformity data	Conformity marking

3.2 Serial Number Label




Fig. 15: Serial number label

A serial number label is affixed inside the housing on the right side of the instrument's back and bears the following data:

Identification data	Model
	REF: Ordering information (material number and revision level)
	SN: Serial number
Address data	Manufacturer's name and address

3.3 Dimensions and Weights



	Dimension	Fluent 480	Fluent 780	Fluent 1080
a	Overall length	1150 mm (45.28 in.)	1650 mm (64.96 in.)	2150 mm (84.65 in.)
b	Footprint depth	780 mm (30.71 in.)		
c	Overall depth	923 mm (36.34 in.)		
d	Overall height on cabinet	1977 mm (77.8 in.)		

Component	Fluent 480	Fluent 780	Fluent 1080
Base unit	120 kg (264.5 lb.)	140 kg (308.6 lb.)	190 kg (418.9 lb.)
Packaging	61 kg (135 lb.)	83 kg (183 lb.)	106 kg (234 lb.)
FCA	10.4 kg (22.9 lb.)		
MCA384	12.6 kg (27.8 lb.)		
384-channel head	7.2 kg (15.9 lb.)		
RGA	10.2 kg (22.4 lb.)		
RGA-Z	10.6 kg (23.4 lb.)		
cXP	1.2 kg (2.6 lb.)		

3.4 Power Supply

NOTICE

Overheating of the Power Supply

The power supply unit can be damaged or destroyed.

- Power supply must not be covered.
- Power supply heat dissipation must be guaranteed.



External devices must not be connected to the power supply. They can lead to a reset or standstill of the Fluent

Tab. 2: Fluent power-input

Supply	Rating
Line voltage (single phase)	100–240 VAC (-15%/+10%)
Input current	9.8 A (at 100 V) – 4 A (at 240 V)
Frequency	50–60 Hz

Tab. 3: Fluent power-output

Supply	Rating
Output voltage	24–28 V factory set: 25.2 V
Continuous power	500 W
Peak power (time limit)	1500 W for 3 seconds
Weight	3.8 kg (8.5 lbs.)

Max. mains supply voltage fluctuation: $\pm 10\%$ of nominal voltage.

Classification with regard to electrical safety according to EN/IEC standards:

Tab. 4: Electrical specifications (safety)

Overvoltage category	II	IEC 60664-1
Pollution degree	2	(EN) IEC 61010-1

3.5 Data and Power Connections




Fig. 16: Data and power connections

A	Fluent instrument	B	Control PC
C	Power supply unit	D	Wall outlet
I1	USB interface	I2	Power cord

The figure shows the components of a sample system with data and power connections. The Fluent instrument parts are shown within the rectangle. The instrument power switch is part of the power supply unit. The power cord is connected to a wall outlet for the mains power supply.

All data traffic to and from the Fluent passes via the USB interface. The USB cable is connected to the PC that controls the instrument.

3.6 Environmental Conditions

CAUTION

Incorrect Pipetting Volumes

Pipetting results can be influenced by operating conditions.

Condensation can influence electronic components.

- If the Fluent is stored or transported at temperatures below room temperature, after installation it will require a few hours for acclimatization.



Operating Conditions

The Fluent is intended for indoor operation and storage only.

Operating temperature	15–32°C (59–90°F)
Operating humidity	30–80% relative (non condensing) at 30°C (86°F)
Operating altitude	max. 2000 m above sea level

Operating conditions for liquid handling and pipetting:

Room temperature	20–25°C (68–77°F)
Operating humidity	30–60% relative (non condensing)
Operating altitude	about 500 m above sea level
Evaporation	An environment with increased airflow (due to laminar flow, air-conditioning or ventilation, etc.), increases the risk of evaporation that can reduce pipetting precision, especially with low volumes or volatile substances. NOTICE! Ensure that validation conditions match the run conditions.

Transport Conditions

Transport temperature	-20 to 60°C (-4 to 140°F)
Transport humidity	20–80% relative (non condensing)

Storage Conditions

Storage temperature	1–60°C (34–140°F)
Storage humidity	5–80% relative (non condensing) at 30°C (86°F) or below

3.7 Emission and Immunity

Noise Emission < 60 dBA (sound pressure), measured at a distance of 1 m from instrument

EMC

The Fluent complies with the emission and immunity requirements described in IEC 61326-1 and IEC 61326-2-6. However, the electromagnetic environment should be evaluated prior to the operation of the Fluent. It is the operator's responsibility to ensure that a compatible electromagnetic environment for the Fluent can be maintained in order that the Fluent will perform as intended.

This equipment is designed for use in a professional healthcare facility environment. It is likely to perform incorrectly if used in a home healthcare facility environment. If it is suspected that performance is affected by electromagnetic interference, correct operation may be restored by increasing the distance between the equipment and the source of the interference.

Do not operate the Fluent in close proximity to sources of strong electromagnetic radiation (e.g., unshielded intentional RF sources), as these can interfere with the proper operation.

4 Description of Function

This chapter explains the basic function of the Fluent, illustrates the structure and provides a functional description of the assemblies.

4.1 Overview

The Fluent is used for pipetting tasks with robotic arms. The robotic arms can aspirate from and dispense to various containers, such as sample tubes or microplates.

The Fluent is available in three different sizes:

- Fluent 480
- Fluent 780
- Fluent 1080




Fig. 17: Instrument overview (instrument may differ from illustration)

- | | | | |
|----------|---------------------|----------|--------------|
| A | Deck | B | Robotic arms |
| C | Options and devices | | |



*An earthquake protection kit for areas prone to earthquakes is also available.
For further information, please refer to section [“Customer Support” \[185\]](#).*

4.2 Deck

Segments

The Fluent deck, which is the sample zone of the instrument, is composed of segments. Deck segments are interchangeable deck components that can have various dimensions and features. **CAUTION! Do not operate the system without deck segments.**

Grid Number

The segment width is expressed in grid numbers. A grid is 25 mm wide and corresponds to the distance between the positioning pins of a segment.

Grid numbers are also used for expressing the location of segments or runners on the deck.

4.2.1 Carriers



Fig. 18: Fluent deck

A Runner

B Segment

Carriers are deck components designed to hold labware or consumables on the deck.

Runners are carriers that slide on and off the grid segments and usually hold sample tubes or reagent troughs.

Segments are static elements locked onto the deck. Some segments have nests (nest segments) that hold labware, such as microplates or deep well plates, or consumables such as DiTi boxes. Some segments have grid pins (grid segments) for loading and unloading runners.

4.2.2 Deck Trays



Fig. 19: Deck tray

Deck trays, which are placed beneath the dynamic deck segments, capture liquid spills that may occur in the manual deck loading area. The system should be operated with as many deck trays as possible installed below the deck to collect all liquid spills. **CAUTION! Do not operate the system without deck trays and deck segments.**

Cut outs in deck trays for tools and instruments are allowed only for the cabinet version.



Fig. 20: Deck trays beneath the deck segments

Deck trays will not be present where the RGA requires access to a device beneath the deck. A set of deck trays is included with the instrument. The deck trays can be washed or replaced as needed. Refer to section [“End of Day”](#) [▶ 108].

4.2.3 Placement on Segment




Fig. 21: Lock pins and positioning pins

A Lock pins

B Positioning pins

Fluent uses pins to position runners, adapters, or options correctly on a segment. The runners are designed to slide onto the pins. Their positioning can then be checked by reading the grid number on the front of the instrument. The lock pins hold the runners in position.

4.2.4 Segment Position




Fig. 22: Side positions and grid positions

A Side positions

B Grid position

The side positions (AB, YZ) can be used to place labware handled by the RGA.



With the FCA or MCA, pipetting is not possible in side positions.

The numbered grid positions (1–n) are accessible to pipetting arms. In multiple arm configurations, however, not all numeric grids are accessible to all pipetting arms. Restrictions may apply depending on the instrument's arm configuration.

4.3 Robotic Arms

Fluent can be equipped with various robotic arms:

- Flexible Channel Arm (FCA)
- Multiple Channel Arm (MCA)
- Robotic Gripper Arm (RGA)

The robotic arms can be equipped with various arm accessories.

4.3.1 Flexible Channel Arm (FCA)



Fig. 23: Flexible Channel Arm



If carry-over is not tolerated, the use of disposable tips with filters is strongly recommended.

The FCA (A) is equipped with pipetting tips and can control liquid handling for up to 8 separate channels.

The FCA configured with DiTi adapters has an optional FCA gripper enabling certain labware moves—refer to [“FCA Gripper”](#) [► 56].

4.3.1.1 FCA with Liquid System (Liquid FCA)

The FCA with a liquid displacement system is filled with system liquid that is supplied by syringe pumps. It is used to pipette liquids with different volume ranges, depending on the tips and the syringe size used. The liquid FCA can be configured with a choice of fixed, washable tips or with disposable tip adapters.



Tecan recommends using deionized water as system liquid.

4.3.1.2 FCA with Air System (Air FCA)

The FCA with an air displacement system is used to pipette liquids by moving a plunger inside the pipetting channel. The Air FCA is configured with disposable tip adapters.

4.3.2 Multiple Channel Arm (MCA)




Fig. 24: Multiple Channel Arm



If carry-over is not tolerated, the use of disposable tips with filters is strongly recommended.

The MCA (A) is a robotic arm with a multiple channel pipetting head. All pipetting head channels aspirate and dispense simultaneously. The pipetting head can exchange head adapters. The different types of head adapters allow various pipetting formats:

- MCA384 with 384 disposable tips
- MCA384 with 96 disposable tips (adapter plate)
- MCA 384 with 384 fixed, washable tips
- MCA 384 with 96 fixed, washable tips

4.3.3 Robotic Gripper Arm (RGA)




Fig. 25: Robotic Gripper Arm

The RGA (A) is a robotic arm with a gripper head and gripper fingers. The RGA transports microplates and other labware between deck positions, peripheral devices and labware storage:

A standard height Robotic Gripper Arm (RGA standard Z) can access objects located on the deck or on the lower deck.

A tall Robotic Gripper Arm (RGA long Z) can access objects located on the deck, on the lower deck and below.

4.3.3.1 Robotic Gripper Head

The RGA can be equipped with two different robotic gripper head options.

The regular gripper head offers a choice of gripper fingers that are manually exchangeable.

The Finger Exchange System (FES) offers automatic finger exchange with a choice of gripper finger sets. Finger sets are mounted on a docking station mounted on a standard nest segment. Finger sets are automatically picked up and placed by the robotic arm. Finger type and finger exchange are monitored. Any or all fingers may be used within a single method.

4.3.4 Arm Accessories




Fig. 26: Arm accessories

- | | | | |
|----------|-----------------|----------|-----------------------------|
| A | Fixed tips | B | Disposable tips |
| C | Gripper fingers | D | Barcode scanner for the RGA |

4.3.4.1 Fixed Tips



If carry-over is not tolerated, the use of disposable tips with filters is strongly recommended.

Washable reusable tips for aspiration and dispense are available for the FCA and MCA.

4.3.4.2 Disposable Tips

Tips are provided in trays or boxes (single or nested) depending on the type. Tips are discarded or re-racked after aspiration. Tips are discarded, with the Disposable Tip Ejection System, into a waste chute mounted on a deck segment.

4.3.4.3 Gripper Fingers

The RGA regular gripper head and the Finger Exchange System (FES) can be equipped with different types of gripper fingers.

Eccentric Gripper Fingers

Eccentric gripper fingers transport microplate-based objects within and beyond the pipetting area. Grasp plate objects from sides. Two variants are available:

- Standard length fingers for loading microplates into hotels and devices.
- Eccentric long fingers for loading deeper devices such as the 4-slot, cell-plate, monitored incubator.

Centric Gripper Fingers

Centric gripper fingers transport microplate-based objects inside and below the pipetting area. Grasp plate objects from above.

Tube Fingers

Tube fingers transport tube-based objects inside and below the pipetting area.

4.3.4.4 Barcode Scanner

The RGA can be equipped with a horizontal scanner for barcodes on microplates and DiTi boxes.



*The laser class safety instructions must be read carefully and must be followed.
Please also refer to the manual provided by the barcode scanner manufacturer.*

4.4 Liquid System (Liquid FCA)



Fig. 27: Liquid system (Liquid FCA)

The liquid system is designed for the efficient washing, inside and outside, of fixed pipetting tips.

4.5 Wash System (MCA)



Fig. 28: Wash system (MCA)

The wash block (A), installed on the MCA segment, washes the tips of the fixed tip adapter after each pipetting cycle.

4.6 Options and Devices



Certain options from Tecan and third-party devices that can be used with Fluent are for research use only (RUO).

In this section research-use-only options and devices are marked with an asterisk (*).

For further information, please refer to section [“Intended Use”](#) [► 11].

Passive Options

- Hotel (plate storage device)
- Cabinet
- Dust cover
- Mix and Pierce

Active Options

- FCA gripper
- HEPA hood
- Fluent Stacker
- MIO2
- Te-Shake
- Te-VacS
- Fluent Carousel
- Wash and refill center (WRC) towers (e.g., MCA wash station)

For further information, please refer to section [“Reference Documents”](#) [► 12].

- Washers based on HydroControl
- Balances based on the MT-SICS level 1 standard
- SiLA-compliant devices*
- Agilent Sealer*

**Barcode
Readers**

Readers

- Inheco ODT
- Inheco Heating Cooling using the MTC/STC controller
- Cytomat 10*, 20*, 200*, and 6000*
- Fluent ID tube barcode scanner
- Barcode readers of the Keyence BL-1300 series
- Tecan Readers controlled by Magellan
- Spark and SparkControl Magellan*
- Ziath 2D flat-bed-reader*



Please also refer to the manuals provided by the option, device or third party device manufacturer. The instructions must be read carefully and must be followed.

4.6.1 Fluent ID Tube Barcode Scanner




Fig. 29: Fluent ID

- | | | | |
|----------|-----------------|----------|-----------------------|
| A | Scanner housing | B | Laser barcode scanner |
| C | Loading area | D | LEDs |
| E | Reflector | | |

The Fluent ID is an optional module that can be incorporated to scan tube barcode labels as tube runners are loaded onto the deck. Each Fluent ID module includes six dedicated grid positions for loading and scanning the barcode labels of up to six runners. The reflector is used to detect empty tube positions in a runner. A graphic interface on the touchscreen monitor provides guidance for the Fluent ID operation.

The laser radiation from the barcode scanner is a low-power, collimated beam in the visible spectrum with the following properties:

- Wavelength: 655 nm
- Pulse duration: 150 μ s
- Maximum power of energy output: 1.0 mW

4.6.1.1 Fluent ID Compatible Tube Runners

The Fluent ID tube runners are each designed for one type of tube:

- Runner with 32 positions for 10 mm diameter tubes
- Runner with 32 positions for 13 mm diameter tubes
- Runner with 26 positions for 16 mm diameter tubes
- Runner with 32 positions for 2ml Eppendorf Safe-Lock tubes



Optional plugs can be used to block two positions of a 26-position runner in order to use it as a 24-position runner, allowing parallel pipetting out of tubes in multiples of eight.

4.6.2 FCA Gripper

Overview




Fig. 30: FCA gripper

The FCA gripper is an option for the FCA configured with DiTi adapters that—in addition to pipetting—allows the FCA to perform some labware moves. The FCA can automatically get and drop the FCA gripper fingers during the run.

FCA Gripper Fingers




Fig. 31: FCA gripper fingers

The FCA gripper fingers need to be replaced after 2 years or 20000 cycles usage (one cycle defined as pickup, use and park). The cycles will be monitored with a counter defined in the Fluent Control software.

FCA Gripper Docking Station Nest




Fig. 32: FCA gripper docking station nest

The FCA gripper docking station nest is used for storing the FCA gripper fingers. It can be mounted like a standard microplate nest on a deck segment.

4.6.3 Mix & Pierce

The Fluent Mix & Pierce Workstation is designed for applications that transfer liquid from and to rubber capped sample tubes without removing the caps by piercing through the rubber caps.

The Fluent Mix & Pierce Workstation is configured with up to 2 liquid FCAs, a deep wash station and up to 4 Tube Rotators, depending on the Fluent base unit size. Tube Rotators can be integrated on any Fluent base unit size and support tube barcode reading, sample mixing, piercing and aliquoting. For more information on supported tube types refer to section [“Tube Rotator Runners” \[► 58\]](#).

The workflow can be divided into the following steps:

1. Barcode scanning while tube loading
2. Mixing of tube contents
3. Piercing and liquid handling with the FCA in the Tube Rotator with Piercing Tips
4. Washing and decontamination of Piercing Tips in the deep wash station and decontamination troughs
5. Repetition of step 2 and continuation

4.6.3.1 Tube Rotator

The main purpose of the Tube Rotator Module is to mix the liquid content of the tubes and to serve as a carrier for piercing and pipetting actions. A single Tube Rotator has a capacity of 5 Tube Rotator Runners with 24 tubes each (i.e., 120 tubes total capacity).

The device contains the following subcomponents:

- An integrated Tube Barcode Scanner in order to scan the sample barcodes during loading
- An optional deep wash station with deep troughs for decontamination of Piercing Tips and an error tube holder. The error tube holder can be used to save samples in case of piercing errors. The wash station is placed next to the rotating drum.
- A rotating drum with tube downholder that accommodates up to five Tube Rotator runners. The drum performs the sample mixing by either 360° rotation or oscillation at different angles and speeds. The downholder (cover) supports the piercing process.
- The Tube Rotator with Tube Rotator Runners support capacitive liquid level detection before and after aspiration as well as after liquid dispense (liquid arrival check) through closed tubes (optional setting).
- The Tube Rotator is installed by the FSE and must not be moved by the key operator or user.

4.6.3.2 Tube Rotator Runners

Tube Rotator runners are designed for the use on the Tube Rotator and to support piercing functionality. There are different tube runners to accommodate the supported tube types for piercing:

- 13x75mm BD Tube Rotator Runner, 24 tube positions
- 13x100mm BD Tube Rotator Runner, 24 tube positions
- 13x75mm Greiner Tube Rotator Runner, 24 tube positions

- 13x100mm Greiner Tube Rotator Runner, 24 tube positions
- 16x100mm Tube Rotator Runner, 24 tube positions

Tab. 5: Tube and runner compatibility

Product Line	Tube			Runner	
	Diameter [mm]	Length [mm]	Labware Definition	Compatible Runner	Bridge Color
Greiner Vacuette	13	100	13x100mm Greiner Vacuette with septum	1x24 13x100mm Greiner Tube Rotator Runner	gray
	13	75	13x75mm Greiner Vacuette with septum	1x24 13x75mm Greiner Tube Rotator Runner	
	16	100	16x100mm Greiner Vacuette with septum	1x24 16x100mm Tube Rotator Runner	black
BD Vacutainer	13	100	13x100mm BD Vacutainer with septum	1x24 13x100mm BD Tube Rotator Runner	white
	13	75	13x75mm BD Vacutainer with septum	1x24 13x75mm BD Tube Rotator Runner	
	16	100	16x100mm BD Vacutainer with septum	1x24 16x100mm Tube Rotator Runner	black

4.6.3.3 Piercing Tip Protection



Fig. 33: Piercing tip protection

The piercing tip protection is a cap that is used for covering the sharp apex of the piercing tips during tip replacement and troubleshooting. It protects the user against injury and the tips against damage.



The piercing tip protection is for one-time use only. After usage all piercing tip protections have to be discarded into the biological waste container.

4.6.3.4 Piercing Tip Removal Tool



Fig. 34: Piercing tip removal tool

The piercing tip removal tool is used for retracting a piercing tip stuck in a tube that cannot be retracted with software commands.

4.6.4 Frida Reader



Fig. 35: Frida Reader

The Frida Reader is intended for automated quantitation and normalization of nucleic acids. Samples to be measured with the Frida Reader need to be cooled at 4°C lest sample evaporation can impair measurement results.

⚠ CAUTION

Vibrations can cause wrong results!

Vibration of the sample drop can cause incorrect measurement results and impair the safety or clinical condition of the patient sample.

- A stable floor is a prerequisite for an appropriate installation location.
- During Frida Reader measurements no internal or external vibration sources are allowed nearby.
- Please avoid sources with the resonance frequency. In particular vibrations around 36 Hz (2160 rpm) and around 42 Hz (2520 rpm) shall be avoided as these are resonance frequencies of a hanging drop.

⚠ CAUTION

Room illumination can cause wrong results!

Room illumination above the module can interfere with the measurement, cause incorrect measurement results, and impair the safety or clinical condition of the patient sample.

- The robotic system has to have a non-transparent cover top, a front and back panel, to prevent ambient light at the measurement position of the Frida Reader.

4.6.5 Phase Separator

The Phase Separator is designed to detect separation phases between liquids of different viscosity. As such, it is independent of the separation phase being visible from outside the labware. The Phase Separator can be used in applications that require the clean transfer of a liquid phase from a source to a destination labware.

The following represents a typical workflow:

- ✓ Centrifugation of the source labware with liquid mixture to generate a distinct phase between liquids.
 - ✓ The liquids must differ in their viscosity to enable phase formation during centrifugation.
1. Load decapped tubes on the Fluent deck and use a barcode scanner (e.g., Fluent ID) for full traceability. Take care not to disturb the layer between the phases/ liquid fractions during loading.
 2. Start the protocol defined for separation of the fractions. The Phase Separator function of the Air FCA will detect the phase between the liquids and start the transfer of fraction of interest to the destination labware.

More than one phase can be extracted from the source labware. The removal of, at least, part of the upper phase is required to prevent spillage of liquid (overflow of labware) on the deck and to prevent potential contamination of the pipetting channel above the disposable tip during phase detection.

For technical specifications refer to the Reference Manual. For details on FluentControl software refer to the Application Software Manual. Refer to ["Reference Documents" \[12\]](#).

The liquid phases must be clearly separated. For the separation of whole blood to plasma and blood cells, the quality and pre-treatment of the samples is essential. Parameters that can have an effect on phase detection in blood samples are sample quality (lipemia, hemolysis), storage time, storage temperature, transportation conditions, centrifugation conditions (time, rcf, temperature, ramp, rotor type), distortion of the phase after centrifugation, etc.

For optimized phase separation results, blood samples should be processed as fast as possible after withdrawal. Sample treatment and storage conditions should follow the specific tube manufacturer's recommendation.

Centrifugation at 2500 rcf for 10 minutes at room temperature with a slowdown ramp lead to a clean phase separation for plasma samples (Tecan internal testing conditions).

To prevent liquid spillage from tubes during aspiration and phase detection, the tubes should not be filled to the rim. Phase detection usually requires a fast downwards movement in the tube combined with a slow aspiration speed, leading to an increase of the liquid level during detection.

For tubes with high fill volume a removal of liquid from the top is recommended before phase detection is started.

5 Control Elements

5.1 Operating Elements



Fig. 36: Operating elements

A Clamp lever

B Touchscreen

Clamp levers lock and unlock the segments.

The touchscreen displays methods and descriptions, allowing the operator to control the instrument.

5.2 User Interface




Fig. 37: User interface of FluentControl


- A** Navigation path **B** Working area
C Display/Option/Action buttons

Through the user interface of FluentControl the operator has access to method runs for operation and system care.

5.2.1 Navigation Path

Use the navigation path to understand and navigate the hierarchy of FluentControl.


Tab. 6: Navigation path buttons

Button	Name	Function
	Home	Press to return to the home page.
	Navigation pane	To display current and previous selections.
	Menu expander	Press the menu expander button to reveal options such as light controls and to switch operators.


5.2.2 Working Area

Access methods and descriptions through the working area of the user interface. Details about the method run status are also displayed here.

Tab. 7: Working area buttons

Button	Name	Function
	Run	Press to start the selected method.
	Add	Press to add more methods to your quick start list.
	Selected method	The method currently selected that will be run when Continue is pressed.
	Available method	A method that can be selected by clicking on it.
	Quick start button	Press to start the selected method immediately.

Tab. 8: Working area display

Display	Display Function
	Displays the status and the remaining time for the method run.
	Description of the currently selected method or additional information on the current action.

5.2.3 Display, Option and Action Buttons


Tab. 9: Display, option and action buttons

Button	Name	Function
	Ok	Press to confirm.
	Cancel	Press to cancel.

Button	Name	Function
	Continue	Press to continue.
	Pause	Press to request a run pause at the end of the current action.
	Stop	Press to halt a run immediately, even in the middle of the current action. If feasible, the system will offer the possibility of restoring or continuing the run.
	Remove	Press to remove the method from the quick start view.
	View mode	Press to toggle between list view and quick start views.
	Sort by	Press to toggle method run display between alphabetical and most recent sorting.

5.2.4 Method Recovery Buttons

Tab. 10: Display, option and action buttons

Button	Name	Function
	Discard	Press to discard a recovered method status.
	—	Press to move to the next screen.
	Recovery Point	Press to return to the previous screen ("Recovery Point").

Button	Name	Function
	Run Recovery	Press to continue the run.

5.2.5 DeckCheck Buttons

Tab. 11: DeckCheck buttons

Button	Name	Function
	Left camera	Displays the camera picture taken from the left camera (Fluent 780/1080 only). An exclamation mark appears on the icon if a layout discrepancy has been seen with this camera.
	Center camera	Displays the camera picture taken from the center overview camera. An exclamation mark appears on the icon if a layout discrepancy has been seen with this camera.
	Right camera	Displays the camera picture taken from the right camera (Fluent 780/1080 only). An exclamation mark appears on the icon if a layout discrepancy has been seen with this camera.
	Pause Alternate	Screen alternates between reference and live pictures: Press this button when either reference or live picture is displayed to hold that picture static.
	Resume Alternate	Picture is static: Press to resume alternating between Reference and Live pictures.
	Check	Activates a recheck of the system—for example when some corrections have been made. Door closure will be prompted. For a 3-arm system, the middle arm must move: If the door is not closed the check will be made, however the middle arm will block one camera.

Button	Name	Function
	Ignore & Continue	Appears only if configured for that command in the method. Allows any highlighted discrepancies to be ignored and the script run will continue.
	Continue	Appears when all discrepancies have been resolved or if the system has not found any discrepancies and the option show always has been selected for the command. This may allow subtle color changes to be seen by eye that the system did not recognize.

Tab. 12: Displays

Display	Description	Function
	Reference picture	The reference picture is stored in the script command displaying the desired deck layout.
	Live picture	The live picture taken by the cameras while running the script.
	Discrepancy (difference to reference picture)	Red squares mark areas where discrepancies have been found between the reference and live pictures. The marked area may include more than one error.

5.3 Error Signals and Instrument Status



Fig. 38: Status lamps

A Power status lamp

B Top status lamp



The status lamps indicate the instrument status by means of different color, steady or flashing lights. The top status lamp is only active when software is running.

Tab. 13: Light signals from status lamps

Signal	Color	Mode	Instrument status
	—	off	The instrument is switched off (disconnected from the power supply).
	white	“heartbeat”	The instrument is switched on (control software connected, modules not initialized yet).
	white (power lamp only)	continuous	Instrument “power on” state (control software is not connected).

Signal	Color	Mode	Instrument status
	color scheme of FluentControl user interface	"heartbeat"	Idle mode All modules are initialized; the instrument is ready to run a method. After about one hour in Idle mode, the instrument will switch to standby mode. Standby mode All axes are braked. The arms are not in ZeroG and cannot be moved manually. To activate the instrument, run a method or request the key operator to select the Move Tool for ZeroG mode.
	yellow	continuous	Teach mode The instrument "learns" positions. In this mode the user can move the robotic arms manually.
	green	continuous	A method (script or process) is running. This is the regular "production" mode.
	red	flashing	Error state The control computer screen or the touchscreen displays an error message.
	color is user configurable	flashing	User prompt System waiting for a user interaction.
	green	flashing	Active stop This is an intentional pause triggered by the runtime controller or by opening a safety panel. The instrument pauses to allow user interaction with the deck. The operator can resume the method.

5.4 Fluent ID Status LEDs



Fig. 39: Fluent ID LEDs

The Fluent ID LEDs signal the following states:

Tab. 14: Fluent ID LEDs

Signal	Color	Mode	Instrument status
	—	off	Fluent ID is idle.
	white	continuous	Fluent ID power on (but not yet initialized).
	blue or custom color	flashing	Ready for runner loading or unloading.
	green	continuous	Barcodes successfully scanned. Runner supervised. Do not unload as this will interrupt the run.
	red	flashing	Error state Error message and required action are displayed on the touchscreen.

6 Operation

6.1 Safety Instructions for This Chapter

CAUTION

Wrong results or contamination of instrument!

Wrong results or contamination of the instrument may occur if the installation qualification and operation qualification have not been performed or if the operating procedures that are given in this manual are not followed.

- Installation qualification and operating qualification records are available and known.
- Methods and processes, including pipetting parameters, must be validated by the key operator.
- Liquid level detection in conjunction with piercing applications for FCA and Air FCA must be validated by the key operator.
- The operator must be trained on the operating procedures, methods and processes.

CAUTION

Biological and chemical contamination of the user!

Damaged FCA gripper fingers can lose plates. Dropped plates can cause contamination by hazardous substances.

- Check the FCA gripper fingers after a crash.

CAUTION

Sharp edges and points!

The Piercing Tips of the Fluent Mix & Pierce Workstation have pointed tips and sharp edges that can cause injuries.

- When loading the instrument, move the FCA to a save position with a software command.
- After an error, cover the piercing tips with piercing tip protections and move the FCA manually to a save position. Refer to section [“Piercing Tip Protection”](#) [▶ 59].

CAUTION

Biological contamination of the system!

In the Fluent Mix & Pierce Workstation, blood can contaminate the caps of the tubes.

- Handle tubes with care.
- Wear protective equipment.

6.2 Operating Modes

Fluent can be run in three different operating modes:

Operator

Routine Operating Mode

- Normal operating mode, where the application or routine system care tasks are run.
- Fluent is monitored by the FluentControl software runtime controller.

Key Operator

Method Definition Mode

- This operating mode is used to perform special tasks such as adjustment to set the method.

FSE

Service Mode

- This operating mode is used to perform special tasks such as tests to ensure the operating readiness.

6.3 Putting into Operation

6.3.1 Switching On the Instrument

To switch on the instrument, proceed as follows:

1. Switch on the power on the power switch (A) on the rear of the external power supply.



When the instrument is powered up, the power lamp will light up blue. Refer to section Error Indication and Instrument Status.

If the status lamp does not light up, start up the PC or contact the key operator.

2. Start up the FluentControl software. Refer to section “Starting FluentControl” [► 74].

6.3.2 Starting FluentControl

- ✓ Operating procedures must be available and known.
 - ✓ Installation qualification and operating qualification records are available and known.
 - ✓ System care has been performed.
 - ✓ Instrument is switched on.
1. Launch the software with **Start > All programs > Tecan > FluentControl**.
*After a few seconds the **Start** screen appears.*


6.3.3 User Login



User Management is only available with the Fluent Gx Assurance Software option.


To log in to FluentControl, proceed as follows:

- ✓ Fluent Gx Assurance Software is installed.
 - ✓ FluentControl has started.
 - ✓ User Management is activated in FluentControl and the process has been defined.
1. Select the assigned user profile (A).



2. Enter password on the keyboard (B).

3. Press **OK (C)**.



Following log-in, the instrument is automatically initialized.

6.3.4 Placing Segments


To place segments, proceed as follows:

- ✓ All segments, carriers, options and devices must be placed in accordance with the selected method.
 - ✓ Segments are cleaned and in perfect condition.
 - ✓ Segments are placed in the corresponding grid position.
1. When installing an MCA Active Carrier, connect the cable.



2. Lower the segment onto the rear of the deck.

3. Align the rear edge to the rear channel cover or the instrument extension.
4. Carefully lower the front part of the segment (A).



5. Turn the clamp lever from left to right, to its closed position. Refer to section [“Checking Segment”](#) [▶ 178].

6.3.5 Removing Segments




Fluent ID segments are not designed to be removed! They are connected directly to the instrument electronics. Connection to electronics can only be made by a qualified FSE.

Please consult section Customer Support.

To remove segments, proceed as follows:


- ✓ All reagents, samples, racks, runners and plates have been removed from the segment.
 - ✓ Nothing is placed on the segment.
1. Turn the clamp lever (A) from right to left to the open position.

The segment is unlocked and the yellow mark on the clamp lever is visible.



2. Push the segment forward by approximately 4 mm.

3. Lift the segment (B) at the front.



Disconnect the cable before removing an MCA Active Carrier.



4. Store the segment in a clean and dry location to avoid any damage.

6.3.6 Loading Standard Runners

NOTICE


Damage due to improper loading or unloading

Damage to runners and pins.


- Align the runner horizontally with the deck.
- Support the front end of the runner with one hand.
- During unloading, ensure that the runner clears all pins before lifting the runner.

To load runners, proceed as follows:

1. Align the runner with the respective grid position (A).



2. Push the runner to the stop position.
3. Ensure that the runner locks the segment securely.
This can be felt in the last couple of millimeters before the runner touches the stop position.



To unload runners, proceed as follows:

1. Pull the runner horizontally at the level of the deck until fully removed from the loading area.
2. Support the front end of the runner with one hand.
3. Ensure that the runner clears all pins before lifting the runner.

6.3.7 Checking the Deck Layout

Ensure that the carriers, labware and devices installed on the deck correspond to the deck layout defined for the method.

NOTICE

Equipment damage!

Incorrect segment and labware positioning on the worktable may cause arms to crash.

- Always ensure that the physical deck configuration and loaded labware matches the FluentControl worktable configuration.
- Always ensure that labware is fitted correctly in the nests. Refer to [“Position Labware”](#) [▶ 153].

NOTICE

Magnetic field creates interference!

A strong magnetic field (north pole up) at the aspiration position may interfere with the tip-presence sensor and may lead to unexpected errors (e.g., **DiTi lost**).

- Ensure that no strong magnet is placed in an SBS-position adjacent to the aspiration position.



Fluent ID segments can only be removed by the FSE because of the connection to the electronic board below the deck.

- ✓ Method must be prepared by the key operator.
 - ✓ The consumables are consistent with the consumables defined in the method.
 - ✓ Fluent ID tube runners must only be loaded after the method has started, when prompted on the touchscreen.
1. Follow the instructions displayed on the touchscreen.

The illustration shows an example of an instruction (A) displayed on the touchscreen:



6.4 Before Starting a Method

The following checklist must to be completed before starting a method.

Tab. 15: Checks before starting a method

Instrument/Component	Task	Reference/Activities
Process validation	Ensure that the method you are selecting has been validated, before starting a production run.	Contact the key operator for further information.
Touchscreen	Follow the instructions on the touchscreen. NOTICE! The instructions provided by the key operator must be strictly observed. If no instructions are displayed, follow the task list below.	—

Instrument/Component	Task	Reference/Activities
Segments, carriers, options and devices	Ensure that all segments, carriers, options and devices are installed and secured. Ensure that only objects intended for use in the method are present on the deck. Ensure that the test run is successfully completed.	If the test run fails, contact the key operator to perform the test run again.
Samples and reagents	Ensure that all samples, reagents and labware are loaded correctly. NOTICE! Barcode scanning only takes place after the method is started. Ensure that the Fluent ID deck is clear of runners before the start of the method. The runners must only be loaded when the prompt is displayed on the touchscreen.	—
Waste tubing (liquid systems only)	Ensure that the waste tubing is routed correctly.	Visually inspect the waste tubing to ensure that they are not kinked or squashed. Replace defective waste tubing. Refer to section “Connecting Wash Station (MCA)” [▶ 152].
Wash system (liquid systems only)	Ensure that the system liquid and waste container are correctly connected.	Refer to section “Checking the Tubing on System Liquid Container and Waste Container” [▶ 84].
Wash system (liquid systems only)	Ensure that the system liquid container is filled to the correct level. Ensure that the waste container is empty.	Refer to section “Connecting the System Liquid Container and Waste Container” [▶ 132].
Wash system (liquid systems only)	Ensure that the correct system liquid is used as defined in the method.	—

Instrument/Component	Task	Reference/Activities
Wash system (MCA wash center only)	Check the liquid level in the wash block.	—
Disposable tip waste and wash station unit	Ensure that the disposable tip waste and the wash station unit are clean.	Refer to section “Cleaning Disposable Tip Waste and Wash Station Unit” [▶ 127].
	Ensure that covers for waste chutes for aerosol containment or for MCA 384 tip guidance are mounted.	—
Disposable tips	Ensure that the correct tips are loaded. Ensure the tip waste is empty.	—
Fixed tips	Ensure that fixed tips are clean and undamaged.	Visually inspect the fixed tips to ensure that they are clean. Visually inspect the fixed tips with a dentist mirror to ensure that the coating is intact.
Deck	Ensure that the carriers, labware and devices installed on the deck correspond to the deck layout defined for the method.	Refer to section “Checking the Deck Layout” [▶ 80].
Labware	Ensure that all labware is positioned securely. If microplates exhibit sideways movement, ensure the labware positioners are correct.	Refer to section “Position Labware” [▶ 153].
Tube Rotator	Ensure that no positioning pins or lock pins are missing on the Tube Rotator.	Refer to section “Replacing Lock Pins and Positioning Pins” [▶ 180]

6.4.1 Checking the Tubing on System Liquid Container and Waste Container

⚠ CAUTION


Sample contamination!

In case of a dual liquid FCA configuration, different types of system liquids can be used for each arm. Connecting the wrong system liquid container to an arm may cause sample contamination.

- Label each system liquid container with the corresponding system liquid name.

✓ The wash system must be correctly installed.

1. Check that the tube (B) is correctly connected to the system liquid container (D).
2. If present, check that the liquid detection sensor (C) is correctly connected to the system liquid container (D).



6.4.2 Checking Waste Container Tubing

1. Check that tube (A) is connected to the waste container (D).



2. Check that the tubes (C) and (D) are correctly connected to the waste container.
3. If present, check that the liquid detection sensor (E) is correctly connected to the waste container.
4. Screw on the lid (F).



6.5 Running a Method

A method is a collection of scripts or processes defined in the FluentControl software. The method can be executed in a run.

The key operator writes a method that can be executed as follows.

NOTICE

Instrument damage!

Instrument damage can result if the deck is not correctly set up or if the software is incorrectly operated or misused.

- Ensure that all safety devices are installed and functional.
- Ensure that the carriers, labware and devices installed on the deck correspond to the deck layout defined for the method.
- Ensure that only objects intended for use in the method are present on the deck.

6.5.1 Starting a Method


- ✓ If applicable (i.e., Fluent Gx Assurance Software is installed and User Management is activated in FluentControl):
Section [“User Login” \[74\]](#) has been performed.
- ✓ Section [“Before Starting a Method” \[81\]](#) has been performed.

1. Select **Method Starter** (A).
Button lights up as soon as it is touched.




2. Select the method (B) to be executed.
Selected method is highlighted.

3. Press **OK**.




4. Press **Run (C)**.




5. Follow the instructions on the touchscreen.
6. If your script includes DeckCheck pay attention to any differences in the actual Live Deck layout compared to the expected Reference Deck layout. Refer to [“DeckCheck Operation” \[97\]](#).
7. Wait for the method run to end.


The screen displays the approximate time (D) when the method run will end.



8. If the **Run finished with error(s)!** message appears, press **Log** to review the errors and warnings.




9. Press **Next** to return to Home screen.



10. Press **Exit**.

The screen below (E) is displayed when the method run ends.



6.5.2 Loading and Unloading Fluent ID Runners



CAUTION

The Fluent is a class 1 laser product pursuant to IEC 60825-1:2014 that emits laser radiation.

Dazzle, flash-blindness and afterimages may be caused by the laser beam.

- Do not stare into the laser beam or into its reflections.

6.5.2.1 Loading Fluent ID Runners

NOTICE

Damage due to improper loading or unloading

Damage to runners and pins.

- Align the runner horizontally with the deck.
- Support the front end of the runner with one hand.
- During unloading, ensure that the runner clears all pins before lifting the runner.

- ✓ Fluent is equipped with a Fluent ID tube barcode scanner.
- ✓ Tubes are loaded in the runners with barcode label facing left.
- ✓ All the tubes in a runner have the same size and shape. For tube runner types refer to section [“Fluent ID Compatible Tube Runners”](#) [▶ 56].

1. Select and start the method using the touchscreen.


*The LEDs start flashing and the message **Please load tubes** is displayed on the touchscreen.*

When using different tube types, ensure that the correct type of runner is used for each specific grid.

2. Support the front end of the runner with one hand.
3. Hold the runner horizontally at the level of the deck.
4. Push the runner to the stop position.



5. Slide the tube runners, one after another, onto the dedicated grids in the Fluent ID loading area.



6. Check that all barcodes were successfully scanned.
The LEDs turn green when the tube runners are in the loaded position and all barcode labels have been scanned successfully.
For Fluent ID LED status description refer to section [“Fluent ID Status LEDs”](#) [71].



7. In the event of a barcode scanning error, unload the runner, correct the problem and load the runner again.
8. Pull the runner horizontally along the deck until it is fully removed.



The Fluent ID reads every code multiple times as it passes the scanner. For small and narrow tubes (i.e., diameter ≤ 10 mm) reduce the speed of manual loading to enable all reads and reduce error reports.



Fig. 40: Barcode reading confirmation displayed on the touchscreen

Tab. 16: GUI meaning (runner)

Square (runner)	Meaning
Green	All tube barcodes in runner read successfully.
White with red outline	Wrong type of runner for this grid position.
Grey with red outline	Missing runner. A runner should be loaded in this grid position.

Tab. 17: GUI meaning (tube position)

Circle (tube position)	Meaning
Green	Barcodes successfully read.
Red	Unreadable barcode
Orange	Duplicate barcode
White with red outline	Missing tube. A tube should be loaded in this position.



When the 2 ml Safe-Lock tube runner is used, it is not possible to differentiate missing tubes from unreadable barcodes. Missing tubes are reported as unreadable barcodes.

6.5.2.2 Unloading Fluent ID Runners

- ✓ The run has finished or a run is in progress and the LEDs are flashing with the message **Please unload tubes** displayed on the touchscreen.
1. Pull the runner horizontally along the deck until it is fully removed.

6.5.3 Loading and Unloading of Tube Rotator Runners

6.5.3.1 Loading of Tube Rotator Runners

CAUTION

Biocontamination of the system and/or user!

Damaged sample tubes can implode leading to spillage of sample on the Tube Rotator.

- Ensure that no damaged tubes are loaded on to the Tube Rotator.


- ✓ Fluent is equipped with a Tube Rotator.
- ✓ Tubes are loaded in the Tube Rotator runners with barcode label facing left.
- ✓ All the tubes in a runner have the same size and shape. For tube runner types refer to section [“Tube Rotator Runners”](#) [▶ 58].

1. Select and start the method using the touchscreen.


*The LEDs start flashing and the message **Please load tubes** is displayed on the touchscreen.*

When using different tube types, please ensure that you chose the correct runner per tube type (either BD or Greiner). Please also ensure that you load tubes of different heights in the corresponding runners: The tubes are always held in position by the runner bridge in the height of their caps. The tube bottoms always have to sit tightly in the tube inserts of the runners.

2. Open the runner locking lever.



3. Support the front end of the runner with one hand.



4. Hold the runner horizontally at the level of the deck.
5. Push the runner to the stop position.
6. Slide the Tube Rotator runners, one after another, onto the dedicated grids on the Tube Rotator.



7. Check that all barcodes were successfully scanned.
The LEDs turn green when the tube runners are in the loaded position and all barcode labels have been scanned successfully.
For Tube Rotator LED status description refer to section [“Fluent ID Status LEDs”](#) [71].
8. In the event of a barcode scanning error, unload the runner, correct the problem and load the runner again.

9. Close the runner locking lever.




Fig. 41: Barcode reading confirmation displayed on the touchscreen

Tab. 18: GUI meaning (runner)

Square (runner)	Meaning
Green	All tube barcodes in runner read successfully.
White with red outline	Wrong type of runner for this grid position.
Grey with red outline	Missing runner. A runner should be loaded in this grid position.

Tab. 19: GUI meaning (tube position)

Circle (tube position)	Meaning
Green	Barcodes successfully read.

Circle (tube position)	Meaning
Red	Unreadable barcode
Orange	Duplicate barcode
White with red outline	Missing tube. A tube should be loaded in this position.

6.5.3.2 Unloading of Tube Rotator Runners



Do not store Tube Rotator Runners that are loaded with tubes outside the liquid handling operating conditions. Refer to section [“Environmental Conditions”](#) [► 44].

- ✓ The run has finished or a run is in progress and the LEDs are flashing with the message **Please unload tubes** displayed on the touchscreen.
- ✓ The Tube Rotator is in horizontal home position.

1. Open the runner locking lever.



2. Pull the runner horizontally along the deck until it is fully removed.



Tube Rotator segments can only be removed by the FSE because of the connection to the electronic board below the deck.

- Method must be prepared by the key operator.
- The consumables are consistent with the consumables defined in the method.
- Tube Rotator Runners must only be loaded after the method has started, when prompted on the touchscreen.

6.5.4 Resetting Errors

If a message is displayed, proceed as follows:

Message

1. Check the display function, button function or error message. Refer to sections [“Working Area” \[▶ 64\]](#) and [“Method Recovery Buttons” \[▶ 66\]](#).
2. Follow the instructions in this manual and on the touchscreen to correct the error.
3. Continue the method run. Refer to section [“Display, Option and Action Buttons” \[▶ 65\]](#).

If the status lamp lights up or changes color, proceed as follows:

Status Lamp

1. Check the instrument status. Refer to section Error Indication and Instrument Status.
2. If the Fluent is equipped with a Fluent ID tube barcode scanner, check the Fluent ID tube barcode scanner LED status. Refer to section [“Fluent ID Status LEDs” \[▶ 71\]](#).
3. Check the display function, button function or error message. Refer to sections [“Working Area” \[▶ 64\]](#) and [“Method Recovery Buttons” \[▶ 66\]](#).
4. Check the troubleshooting table. Refer to section [“Troubleshooting Tables” \[▶ 141\]](#).
5. If the problem cannot be solved, please consult [“Customer Support” \[▶ 185\]](#).

6.6 DeckCheck Operation

If your script includes use of the DeckCheck the DeckCheck camera system will take pictures of the worktable after loading and compare the actual Live layout to the Reference layout.

The DeckCheck takes approximately 20 seconds for a 3-arm/ 3-camera system and approximately 12 seconds for 1- or 2-arm/ single camera system to take pictures of the deck and display the comparison of the live and reference layouts (assuming that the PC configuration is appropriate—refer to FluentControl Application Software Manual).


Note that for the first use after Instrument power only, the DeckCheck command will need longer to display first result—this may take a few minutes.

During this time the rear LED will be switched on.

On 3-arm Fluent systems the middle arm needs to move between left and right positions (on 1- or 2-arm systems, the left and right arms will be positioned on the far left and far right sides respectively.) For this arm move the front door must be closed. If the image is taken while the door is open, one camera will generally be blocked by the middle arm or by any arm that has been moved manually.

During the DeckCheck process the touchscreen displays shapes moving across the screen and **Taking Images** followed by **Checking**. After 12–20 seconds depending on the instrument size and configuration the deck images will be displayed in alternating mode. The image displayed will be the first camera with a noted discrepancy starting from the left.

DeckCheck screen displaying a discrepancy relative to the Reference picture. Here, the Reference picture shows that a plate should be present and the center camera has detected the discrepancy.



If you are offered the **Ignore & Continue** bottom, the run will continue with the prevailing worktable. Select the **Ignore & Continue** button if you are certain there are no more differences to the required worktable and before closing the door. Select **Check** if you wish to take new images of the Deck—note that if the door is not closed on a 3-arm systems the image will be taken but the middle arm will obstruct the camera. A check is otherwise automatically executed when the door is closed.



Note that if the script includes the option **show always**, the screen above will be displayed and no discrepancies are highlighted. However, the reference and live pictures will alternate and it may be that there are small differences that are not captured by the system but are easily seen by eye—for example some color differences, single missing tubes/tips or small lateral shifts. Refer to the limits listed below.

If discrepancies are detected these will be highlighted.

To correct differences:

1. Open the door and replace or correct the position of the highlighted items.
2. The DeckCheck will work continuously to compare the corrected Live situation with the Reference layout.
3. Use the DeckCheck buttons to look at differences captured by each of the cameras or, to pause the view and hold the Reference picture as needed.
When no further differences are detected the green continue button will appear.
4. Select **Continue** to proceed with the method.



*If any remaining differences are in fact judged acceptable (e.g., total number of tips may be variable at the start of the method or liquid levels vary significantly at the beginning of the run) you may select **Ignore & Continue** if offered in the script by your key operator.*

Some layout differences may not be highlighted by the DeckCheck—e.g., the following colored FCA tip trays:

Difference between MCA head adapter types:

- Yellow/ orange
- White/ orange
- Grey/ all colors

MCA 384 different tip types

MCA 96 different tip types

Missing tubes on partially loaded tube runners

Trough 300 SBS

Microplates rotated through 180 degrees

Microplate well shape (e.g., round versus flat bottomed or PCR well)

Plates in peripheral hotels 10 ml/ 25 ml troughs as insert

Some transparent lids

Many of these differences are however clearly visible in the switch between live and reference layouts.

6.7 Method Recovery

FluentControl offers the option to recover from errors—e.g.:

Previous method run was aborted or had a fatal error: The method recovery option offers the possibility of continuing from the point at which failure occurred in the previous run.



After a method was aborted or had a fatal error, daily maintenance shall be executed. Refer to “Daily System Care” [107].

6.7.1 Switching to Method Recovery Mode

- ✓ Key operator has enabled the method recovery option in FluentControl.
- ✓ Previous method run was aborted.

1. Select **Method Recovery** (A).




6.7.2 Recovering a Method Run

- ✓ Section “Switching to Method Recovery Mode” [▶ 100] has been performed.


1. Select continue to the next screen (A).

The screen displays the last script line executed (C) and the script line where failure occurred—the recovery point (B).



2. Select required buttons described in section “Method Recovery Buttons” [▶ 66].

3. Ensure that the physical deck layout of the Fluent matches the worktable layout (D) displayed on the touchscreen.
 4. Select **Run Recovery** (E).
- The system will start.*



6.8 Switching Off the Instrument

If no method is running, the instrument switches to standby mode. The instrument does not need to be switched off at the mains.

To switch off the instrument, proceed as follows:

1. Stop any method and select the standby mode on the touchscreen.
2. Ensure the instrument is in standby mode. Refer to section Error Indication and Instrument Status.
3. Place the robotic arms in a move free area.

4. Switch off the power using the power switch (A) on the rear of the external power supply.



7 System Care

This chapter gives instructions on all system care tasks to be performed to maintain the Fluent in good working condition.



Only operate the Fluent when it is in good working condition. Strictly observe the system care instructions as described in this manual. To ensure optimum performance and reliability, perform the maintenance and cleaning tasks regularly.

In the event of any problems and for inquiries please consult section “[Customer Support](#)” [▶ 185].

7.1 Decontamination

Decontamination, according to standard laboratory regulations, is required under the circumstances listed in section “[Decontamination Declaration](#)” [▶ 38].

WARNING

Contamination!

Substance residues on the Fluent can cause personal injury and affect the integrity of the process.

- Decontaminate the Fluent and its parts and accessories before any interaction.

The decontamination method must be defined by the key operator based on the type of contaminant and degree of contamination. Guidance on the selection of decontamination agents and application modes is provided in this chapter.



For information on Hydrogen Peroxide Vapor treatment refer to the Reference Manual. Refer to “[Reference Documents](#)” [▶ 12].

CAUTION

Incorrect measurement results of the Frida Reader!

If no insert is mounted, the Frida Reader may deliver incorrect measurement results.

- Use the red blind plug, if the insert is taken out (e.g., for cleaning).

7.2 Cleaning Agents

7.2.1 Cleaning Agents Specifications

Special cleaning agents are required for system care. All the recommended cleaning agents have been carefully selected and tested.

NOTICE

Reduced effectiveness and chemical compatibility!

There is no guarantee for the effectiveness of cleaning agents and chemical compatibility if other cleaning agents than those recommended by Tecan are used.

- Only use cleaning agents recommended by Tecan.
- Cleaning agents are defined for each specific use in the system care tables. Do not use cleaning agents if not specified for use for a specific task.

The following table specifies the cleaning agents referred to in this manual:

Tab. 20: Cleaning agents

Agent	Specification
DI water	Distilled or deionized water
Alcohol	70% ethanol, 100% isopropanol (2-propanol)
Weak detergent	Liqui-Nox
Surface active agent	Contrad 70, Contrad 90 / Contrad 2000, Decon 90
Disinfectant	Bacillol plus, SporGon
Surface disinfectant (for nucleic acid contamination)	DNAzap
Weak acid	sulfuric acid 0.3M, 10% acetic acid, 30-40% formic acid
Base	sodium hydroxide 0.1M
Bleach	2% sodium hypochlorite
System liquid	As defined in the method. Note that Aqueous solutions with salt content should be flushed out during system inactivity—e.g., overnight or weekends. See System Care “End of Day” [► 108].

7.2.2 Commercial Cleaning Agents

All instructions—given by the manufacturer of the cleaning agents or provided in this manual—for handling the cleaning agents must be carefully read and followed.

The table below lists a number of commercially available cleaning agents and disinfectants:

Tab. 21: Commercial cleaning agents

Cleaning Agent	Agent category	Manufacturer
DNAzap	Surface disinfectant (for surfaces contaminated with nucleic acids)	Ambion www.ambion.com
Decon, Contrad	Surface active agent	Decon Laboratories www.deconlabs.com
SporGon	Disinfectant	Decon Laboratories www.deconlabs.com
Bacillol Plus	Disinfectant	www.bode-chemie.com
Liqui-Nox	Weak detergent	Alconox www.alconox.com


7.3 System Care Mode

The key operator defines the system care methods required according to the system care tables provided in “[System Care Tables](#)” [106]. The **System Care** mode, accessed on the touchscreen, provides guidance for system care tasks.


7.3.1 Switching to System Care Mode

✓ System care methods must be available.

1. Select **System Care** (A).



2. Select the task to be executed.



3. Press **Run** to initiate the system care method.
4. Carry out the system care tasks.

7.3.2 Resetting Errors

If a message is displayed, proceed as follows:

Message

1. Check the display function, button function or error message. Refer to sections [“Working Area” \[▶ 64\]](#) and [“Method Recovery Buttons” \[▶ 66\]](#).
2. Follow the instructions in this manual and on the touchscreen to correct the error.
3. Continue the method run. Refer to section [“Display, Option and Action Buttons” \[▶ 65\]](#).

If the status lamp lights up or changes color, proceed as follows:

Status Lamp

1. Check the instrument status. Refer to section Error Indication and Instrument Status.
2. If the Fluent is equipped with a Fluent ID tube barcode scanner, check the Fluent ID tube barcode scanner LED status. Refer to section [“Fluent ID Status LEDs” \[▶ 71\]](#).
3. Check the display function, button function or error message. Refer to sections [“Working Area” \[▶ 64\]](#) and [“Method Recovery Buttons” \[▶ 66\]](#).
4. Check the troubleshooting table. Refer to section [“Troubleshooting Tables” \[▶ 141\]](#).
5. If the problem cannot be solved, please consult [“Customer Support” \[▶ 185\]](#).

7.4 System Care Tables

To ensure optimum performance and reliability, perform the maintenance and cleaning tasks as recommended.



The tasks in the system care tables can only be carried out in the System Care mode. Refer to section [“System Care Mode” \[▶ 105\]](#).

The system care tasks must be performed at regular intervals—namely, daily, weekly, and monthly system care.

7.4.1 Daily System Care

7.4.1.1 Beginning of Day

Run the **DailySystemCare** method, if made available by the key operator; or perform each individual task, applicable to your Fluent arm configuration, listed in the table below in chronological order.

Tab. 22: Beginning of day system care table

Instrument/ Component	System Care Task	Cleaning Agent/ Disposable Product/ Device	Reference/ System Care Activities
Piercing tips	Visually inspect the piercing tips for deposits. Clean if necessary. Inspect that tips are not bent.	70% ethanol or 2% bleach and lint-free cloth	Refer to section “Cleaning Piercing Tips” [▶ 118].
Disposable tip cones and fixed tips	Inspect for damage and deposits	–	This task is included in the Daily System Care method. NOTICE! Disposable tips are not intended for reuse.
Fixed tips	Clean. Inspect that tips are not bent. Visually inspect with a dentist mirror to ensure that the coating is intact.	70% ethanol or 100% isopropanol and lint-free cloth	Refer to section “Cleaning Fixed Tips” [▶ 118].
System liquid container (Liquid FCA and MCA with fixed tips)	Ensure it is clean and full with no visible bubbles Ensure that the tubing to container connectors are joined properly	–	This task is included in the Daily System Care method.
Liquid waste container (Liquid FCA and MCA with fixed tips)	Ensure that it is empty Ensure that the tubing to container connectors are re-joined properly	–	This task is included in the Daily System Care method.
Disposable tip waste bag	Ensure that it is empty	–	Refer to section “Changing Disposable Tip Waste Bag” [▶ 130]. This task is included in the Daily System Care method.

Instrument/ Component	System Care Task	Cleaning Agent/ Disposable Product/ Device	Reference/ System Care Activities
Liquid system (Liquid FCA)	Ensure that it is clean	System liquid, alcohol, DI water	This task is included in the Daily System Care method or can be run separately as Liquid FCA Routine Flush Maintenance method. Refer to section “Cleaning Liquid Path” [▶ 131].
Wash system (MCA with fixed tips)	Flush/prime	Wash liquid	Run the direct command Prime Wash Station (MCA384) .
Liquid system (Liquid FCA)	Visually check for absence of droplets on the tips or DiTi cone after flushing	–	This task is included in the Daily System Care method.
Gripper fingers	Check that fingers are straight and level Inspect for damage and misalignment	–	In case of misalignment refer to section “Robotic Gripper Arm (RGA) Troubleshooting” [▶ 148]. Deformation or damage. Please consult section “Customer Support” [▶ 185].
FCA gripper fingers	Inspect for damage	–	Replace if damaged. For ordering information refer to the Reference Manual. Refer to section “Reference Documents” [▶ 12].
Frida Reader	Remove the blind plug and fit the insert into the Frida Reader	–	Refer to section “Frida Reader” [▶ 139].

7.4.1.2 End of Day

The following table lists the daily system care tasks at the end of the day in chronological order:

Tab. 23: End of day system care table

Instrument/ Component	System Care Task	Cleaning Agent/ Disposable Product/ Device	Reference/ System Care Activities
Deck trays	Check for spills and clean or replace as needed.	DI water, alcohol, weak detergent, disinfectant, base, bleach, DNAzap	Refer to section “ Cleaning Deck Trays ” [▶ 123].
Segments Fluent ID housing	Clean	DI water, alcohol, weak detergent, disinfectant, base, bleach, DNAzap NOTICE! The scanner window requires different cleaning agents than the segment itself. Refer to “ Weekly System Care ” [▶ 111].	Refer to section “ Cleaning Runners and Segments ” [▶ 123]. WARNING! Do not stare into the laser beam.
Reflector foil (Fluent ID, Tube Rotator)	Clean and inspect for damage	Alcohol NOTICE! The reflector foil requires different cleaning agents than the segment itself.	Damage. Refer to section “ Replacing Fluent ID Reflector Foil ” [▶ 125].
Runners	Clean	DI water, alcohol, weak detergent, disinfectant, surface active agent, weak acid, base, bleach, DNAzap	Refer to section “ Cleaning Runners and Segments ” [▶ 123].
Fixed tip block (MCA)	Park, clean, check and cover	Alcohol	Refer to section “ Cleaning Fixed Tip Block (MCA) ” [▶ 138]. CAUTION! Tecan recommends you to park the tip block and to remove it from the instrument at the end of every work shift. Always clean and properly store the tip block. Check for tip displacement.
Fixed tips	Clean	Alcohol, bleach, lint-free cloth	Refer to section “ Cleaning Fixed Tips ” [▶ 118].
Disposable tip cones	Clean	Alcohol, lint-free cloth	Refer to section “ Cleaning Disposable Tip Cone ” [▶ 117].

Instrument/ Component	System Care Task	Cleaning Agent/ Disposable Product/ Device	Reference/ System Care Activities
Wash and waste station (Liquid FCA)	Clean	DI water, alcohol, weak detergent, disinfectant	Refer to section “ Cleaning Disposable Tip Waste and Wash Station Unit ” [▶ 127].
Disposable tip waste slide and waste covers	Clean	DI water, alcohol, weak detergent, disinfectant	Refer to section “ Cleaning Disposable Tip Waste Slide ” [▶ 128].
Liquid system (Liquid FCA)	Flush	System Liquid NOTICE! If the liquid system has a high salt content, flush with deionized water.	Run the Liquid FCA Routine Flush Maintenance method.
MCA wash station (MCA with fixed tips)	Clean	Alcohol, bleach, Decon 90, Contrad NOTICE! Surface active agents, such as Decon or Contrad can affect the process. Therefore, if such agents are used, carefully validate the process.	Run the direct command Prime Wash Station (MCA384) .
Disposable tip waste bag	Change	Recommended bag specifications: W x L: 300 mm x 600 mm; Thickness: 0.5 mm Material: Polypropylene, Polyethylene or copolymer (autoclavable) NOTICE! The waste bag used must comply with the local safety guidelines.	Refer to section “ Cleaning Disposable Tip Waste Slide ” [▶ 128].
System liquid container (Liquid FCA)	Ensure that it is clean	System liquid	Refer to section “ Connecting the System Liquid Container and Waste Container ” [▶ 132].

Instrument/ Component	System Care Task	Cleaning Agent/ Disposable Product/ Device	Reference/ System Care Activities
Waste container (Liquid FCA and MCA with fixed tips)	Empty and clean	DI water, alcohol, weak detergent, surface ac- tive agent, disinfectant, base, bleach	Depending on your local lab- oratory rules/regulations clean daily or weekly. Refer to section “Connecting the System Liquid Container and Waste Container” [▶ 132].
Safety panel	Clean	DI water, alcohol, weak detergent	Refer to section “Cleaning Safety Panels” [▶ 127].
Tube Rotator	Clean surfaces, down- holder and wash sta- tion	Linth-free cloths with 2% bleach, 70% ethanol or 100% iso- propanol	Refer to section “Cleaning the Tube Rotator” [▶ 118].
Piercing tips	Visually inspect the piercing tips for de- posits. Clean if neces- sary.	70% ethanol or 2% bleach and lint-free cloth	Refer to section “Cleaning Piercing Tips” [▶ 118].
Frida Reader	Remove the insert and fit the blind plug into the Frida Reader	–	Refer to section “Frida Reader” [▶ 139].

7.4.2 Weekly System Care

Weekly system care should be performed on the last working day of each week.

Run the **WeeklySystemCare** method, if made available by the key operator; or perform, in addition to the daily tasks, each individual task, applicable to your Fluent arm configuration, listed in the table below in chronological order.

Tab. 24: Weekly system care table

Instrument/ Component	System Care Task	Cleaning Agent/ Disposable Product/ Device	Reference/ System Care Activities
Wash system (MCA with fixed tips)	Visually check for dirt in the tubing and filter	Filter	Dirty filter. Please consult section “Customer Support” [▶ 185].
Pipetting head (MCA)	Perform a leakage test	–	Run the MCA384 leakage and zero dispense method.
Plate adapter (MCA)	Clean	Alcohol Compressed air (for drying)	Refer to section “Cleaning Plate Adapter (MCA)” [▶ 138].

Instrument/ Component	System Care Task	Cleaning Agent/ Disposable Product/ Device	Reference/ System Care Activities
Liquid system (Liquid FCA)	Clean	Depending on the liquid handled by Fluent Decon, Contrad, base, weak acid, disinfectant Followed by water, alcohol and system liquid flushes	Refer to section “Cleaning Liquid Path” [► 131].
Liquid FCA	Check the correct tightness of the syringes at the valve interface and the correct tightness of the syringe plunger at the plunger lock screw.	—	Refer to section “Checking Tightness of Syringes” [► 133]
Liquid FCA	Perform a leakage test (Liquid FCA)	—	Run the FCA Leakage Method .
Air FCA	Perform a leakage test (Air FCA)	—	Run the Air FCA Leakage Method .
System liquid container	Clean	DI water, alcohol, weak detergent, surface active agent, disinfectant, base, bleach	Refer to section “Cleaning the System Liquid Container and Waste Container” [► 133].
Wash station (Liquid FCA)	Clean	Detergent or antiseptic solution	—
RGA gripper finger pads	Remove particles and residues from gripper finger pads	Lint free cloth with alcohol	Wiping with cleaning agent
Docking station and Gripper fingers (attachment interface)	Remove particles and residues from gripper finger attachment interface (PCBA, magnet and conus)	Lint free cloth with alcohol	Wiping with cleaning agent
Stand-alone barcode scanner window	Clean	Weak detergent	WARNING! Do not stare into the laser beam. Refer to the barcode scanner manufacturer's manual. Refer to section “Laser Radiation Instrument” [► 36].

Instrument/ Component	System Care Task	Cleaning Agent/ Disposable Product/ Device	Reference/ System Care Activities
Fluent ID and Tube Rotator scanner window	Check for dirt and damage Clean if necessary	Weak detergent DI water for rinsing	WARNING! Do not stare into the laser beam. Clean and rinse using a soft cloth.
Fluent ID and Tube Rotator reflector	Check for dirt and damage Clean if necessary	Weak detergent DI water for rinsing	WARNING! Do not stare into the laser beam. Clean and rinse using a soft cloth.
FCA gripper	Clean	Alcohol	–
DiTi cones	Check tightness of DiTi cone	–	Refer to section “Tightening a DiTi Cone” [138]
Tube Rotator	Check presence and tightness of lock and positioning pins. Tighten or replace the pins if necessary	–	Refer to section “Replacing Lock Pins and Positioning Pins” [180] .

7.4.3 Monthly System Care

The following table lists the monthly system care tasks in chronological order:

Tab. 25: Monthly system care table

Instrument/ Component	System Care Task	Cleaning Agent/ Disposable Product/ Device	Reference/ System Care Activities
Software	Restart the computer	–	Switch computer off. Wait 10 seconds. Switch computer on again.
Arm guide	Clean	Cotton swab or a lint-free cloth on a screwdriver	Refer to section “Cleaning Arm Guide” [138] .

7.4.4 Periodic System Care



The intervals at which these tasks must be performed should be determined by the key operator.

The following table lists the system care tasks in chronological order:

Tab. 26: Periodic system care table

Instrument/ Component	System Care Task	Cleaning Agent/ Disposable Product/ Device	Reference/ System Care Activities
Cone-sleeve connection	Remove particles Clean surfaces	Alcohol, lint-free cloth	–
UVC light	Check for fingerprints. Clean if necessary.	Alcohol, lint-free cloth	

7.4.5 Yearly System Care

The yearly system care helps to maintain the accuracy and precision and to minimize instrument downtime. It also helps to prolong the life-span of the Fluent.

Please contact the local Tecan service organization to schedule the yearly system care appointment. Please consult section [“Customer Support”](#) [▶ 185].

7.4.6 Biennial System Care

The following preventive maintenance tasks must be performed every 2 years:

Tab. 27: Biennial system care

Component	Task	Reference
FCA gripper	Replace FCA gripper fingers. Reset counter in FluentControl.	For ordering information refer to the Reference Manual. Refer to “Reference Documents” [▶ 12].

7.5 System Care Activities

To perform the system care activities described below, proceed as follows:

- Switch to System Care mode. Refer to section [“System Care Mode”](#) [▶ 105].
- Follow the instructions as described below.

7.5.1 Moving the Instrument on a Cabinet within the Laboratory

CAUTION

Damage to the cabinet!


Cabinet shelves may have been removed, for example, for centrifuge installation. Moving the instrument placed on a cabinet without installed cabinet shelves may damage the cabinet and cause injuries.

- Before moving the instrument, install the cabinet shelves.


To move the instrument on a cabinet within a room, proceed as follows:

1. Ensure that the cabinet is safely parked and secured against rolling away.

2. Ensure that the cabinet shelves (A) are installed.



3. Turn the nut on the cabinet feet (B), using an open-end spanner.
4. Turn the red screw (D) on the cabinet feet (C) until the lock is released and the wheels are in the moving position.




5. Move the instrument on the cabinet to the new location.
6. Ensure that the cabinet is safely parked and secured against rolling away.

7.5.1.1 Leveling the Instrument

To level the instrument, proceed as follows:

1. Using an open-end wrench, lower all adjustable feet until the cabinet wheels can be rotated by hand.




2. Loosen the locknut (A) on the corresponding foot.



3. Place the reference segment according to the grid positions listed below.
Instrument size 480: Grid position left side 1 and grid position right side 21.
Instrument size 780: Grid position left side 1 and grid position right side 41.

Instrument size 1080: Grid position left side 1 and grid position right side 59.



4. Use the spirit level to ensure that the instrument is horizontally and vertically aligned.
5. Adjust the cabinet level as required (clockwise to raise, counterclockwise to lower).



6. After leveling the instrument, retighten the locknuts on the cabinet feet.
7. Ensure that the cabinet is safely parked and secured against rolling away.

7.5.2 Cleaning Disposable Tip Cone

To clean the disposable tip cone, proceed as follows:

1. Clean the disposable tip cones with alcohol, using a lint-free cloth.
2. Check the disposable tip cones and the protruding tip during system care.

For liquid FCA: Ensure that the tubing extension protruding outside of the cone is not damaged.

3. Ensure that the tubing extensions are clean and free of deposits.

7.5.3 Cleaning Fixed Tips

CAUTION

Risk of injury from fixed tips during cleaning

Pipetting fixed tips can cause injuries.

- Avoid contact with the pipetting tips and contact with aerosols when accessing the worktable, by wearing adequate protective clothing.

To clean the fixed tips, proceed as follows:

1. Clean the fixed tips with alcohol, using a lint-free cloth.
2. Ensure that the fixed tips are clean and free of deposits.

7.5.4 Cleaning Piercing Tips

To clean the piercing tips run the **Piercing Tip Cleaning Maintenance** method. This method needs to be adjusted according to your worktable setup.

The script includes the following steps:

1. Prepare the worktable (i.e., labware and hardware).
2. Pierce to Z-start of 8 empty capped tubes on a Tube Rotator or on a tube downholder carrier.
3. Manually clean the accessible part of the piercing tips with 70% ethanol or 2% bleach using a lint-free cloth. Avoid contact with the sharp apex of the piercing tips.
4. Perform wash commands after manual cleaning.


7.5.5 Cleaning the Tube Rotator

General Cleaning Procedure

1. To clean any parts of Tube Rotator use lint-free cloths and soak them with one of the following cleaning liquids: 2% bleach, 70% ethanol, 100% isopropanol
2. Wipe the parts with the soaked cloths to clean and disinfect.
Use cotton swabs to clean areas that cannot be reached with a lint-free cloth.
3. Wipe off the cleaning liquids with cloths soaked with water within 5 minutes after applying the cleaning liquids.

Removal and Cleaning of Tube Downholder Plate

1. To release the downholder plate (A) hold it with one hand and pull the downholder locking pin (B) with the other hand.



2. Remove the downholder plate from the tube rotator.




3. Clean the downholder plate according to the general instructions above or alternatively the tube downholder can be incubated in a bath with 2% bleach for maximal 2 hours.

Cleaning of the Tube Rotator Surfaces

1. Clean the accessible surfaces of the Tube Rotator according to the general instructions above.

2. To change the drum position manually, hold the drum with one hand and press the solenoid release button.



3. Rotate the drum manually and release the solenoid release button.
4. Rotate the drum until it is locked by the solenoid.
5. Clean the surfaces that were not accessible before according to the general instructions above.

Mounting of the Downholder Plate

1. Place the downholder plate on top of the tube rotator drum.
2. Press the downholder plate towards the bottom of the instrument with one hand and push the black slider to the back in order to lock the downholder plate into position.




7.5.6 Cleaning the Tube Rotator Wash Station

General Cleaning Procedure

- ✓ The wash station can be cleaned on the worktable or it can be unmounted for cleaning.
 - ✓ Use a bottle brush instead of cloths for better access.
1. To clean any parts of wash station use lint-free cloths and soak them with one of the following cleaning liquids: 2% bleach, 70% ethanol, 100% isopropanol
 2. Wipe the parts with the soaked cloths to clean and disinfect.
 3. Wipe off the cleaning liquids with cloths soaked with water within 5 minutes after applying the cleaning liquids.

Unmounting of Tube Rotator Wash Station

1. Press the wash station release lever towards the barcode scanner housing and lift the wash station with the other hand.



2. Disconnect the waste tubing and place the connectors into the waste tubing holders.




Mounting of Tube Rotator Wash Station

1. Connect the waste tubing connectors.



2. Mount the wash station on to the guidance shafts (A) and press it onto the base plate.
Verify that the release lever snaps back into place and holds the wash station in position.



7.5.7 Cleaning Runners and Segments

To clean the runners and segments, proceed as follows:

NOTICE

Malfunction of liquid detection (cLLD)!

Possible malfunction of the liquid detection (cLLD) due to compromised contact between runner and deck segment.

Always ensure that the runners and segments are clean and dry.

1. Remove the runners from the instrument deck.
Segments and nests are getting cleaned in place.
2. Wipe the surface of the runners, segments and nests with the cleaning agent.
Rinse the runners, segments and nests with DI water.
3. Place the runners back on the instrument deck.

7.5.8 Cleaning Deck Trays

To clean the deck trays, proceed as follows:

- ✓ Segments above the deck tray are removed. Refer to [“Removing Segments”](#) [76].
 - ✓ If deck segments, such as Fluent ID or MCA Active Carrier cannot be removed, slide deck trays to an open deck position.
1. Remove the deck trays from the instrument.

2. Empty the trays by removing the liquid according to laboratory handling protocol for that liquid.
3. If the deck trays are damaged or lost, they must be replaced.
4. Wipe surface of the deck trays with the cleaning agent.
5. Place the deck trays back in the instrument.

Orient the deck trays as depicted below.

Adjacent deck trays must interlock.




Fig. 42: Incorrect deck tray placement




Fig. 43: Interlocking deck trays





Fig. 44: Correct deck tray placement

7.5.9 Replacing Fluent ID Reflector Foil

✓ Self-adhesive reflector foil

1. Heat up the reflector foil. Use a heat gun.
2. Remove the reflector foil.



3. Remove any residues with alcohol.

4. Apply the new self-adhesive reflector foil on the upper end of the reflector.



7.5.10 Applying Fluent ID Reflector Foil on DiTi Waste Slide

- ✓ Self-adhesive reflector foil

1. Apply tape on the DiTi waste slide according to the illustration below.



2. Apply the new self-adhesive reflector foil on the DiTi waste slide according to the illustration below.

The laser beam must be in the center of the reflector foil.



3. Remove the tape from the DiTi waste slide.

7.5.11 Cleaning Safety Panels

To clean the safety panels, proceed as follows:

1. Wipe the inner and outer surface of the safety panels with cleaning agent.

7.5.12 Cleaning Disposable Tip Waste and Wash Station Unit

To clean the disposable tip waste and the wash station unit, proceed as follows:

1. Push the quick-release fastener button (B).
2. Slide the wash station backwards.




Fig. 45: Fastener for bag housing removal

3. Remove the wash station from the disposable tip waste and the wash station unit.




Fig. 46: Remove the wash station

4. Wipe the surface of the wash station with a cleaning agent and remove any spilled reagents.
5. Push the quick-release fastener button (B).
6. Place the wash station (A) in position.
7. Push the wash station forward.




Fig. 47: Refit the wash station

7.5.13 Cleaning Disposable Tip Waste Slide

To clean the disposable tip waste slide, proceed as follows:

- ✓ Front safety panel is open.
1. Remove the cover (A) from the disposable tip waste slide.




Fig. 48: Removing the cover from disposable tip waste slide

2. Remove the disposable tip waste slide (B) from the holder.




Fig. 49: Removing the disposable tip waste slide

3. Hold a tissue under the disposable tip waste slide bottom opening (C).
Prevent dripping of contaminated substances.




Fig. 50: Handling of disposable tip waste slide

4. Place the disposable tip waste slide and the cover in a basin filled with cleaning agent.
5. Leave to soak for 30 minutes to 4 hours.
6. Take the disposable tip waste slide and cover out of the basin and place them on a clean dry towel.

7. Leave to dry.
8. Reinstall the disposable tip waste slide (B) on the holder.




Fig. 51: Reinstall disposable tip waste slide inset

9. Ensure that the positioning pin is correctly inserted in the slot (D).
10. Place the cover (A) on top of the waste slide.




Fig. 52: Positioning pin and cover

7.5.14 Changing Disposable Tip Waste Bag

To change the disposable tip waste bag, proceed as follows:

1. Lift the fastener (A) and slide the bag housing forward.




Fig. 53: Fastener for bag housing removal

2. Remove the disposable bag housing (A).
3. Remove the disposable tip waste bag (B).
4. Dispose of the disposable tip waste bag according to your laboratory guidelines.
5. Insert the new disposable waste bag (B) into the empty bag housing (B).




Fig. 54: Bag housing and disposable tip waste bag

6. Slide the bag housing into position and close with fastener (A).




Fig. 55: Closing fastener

7.5.15 Cleaning Liquid Path

- To remove protein residues from the inside of fixed tips, use weak acid followed by a base cleaning agent.
- To remove nucleic acid residues from the inside of fixed tips, use a base cleaning agent.
- Cleaning agents, such as Decon/Contrad can affect the process. Therefore, if such agents are used, carefully validate the process.
- Isopropanol is a highly effective disinfectant. It evaporates quickly, leaving surfaces ready to use.
- Only use allowed cleaning agents. Do not use bleach solutions to flush the entire liquid system.

To clean the liquid path, proceed as follows:

1. Detach the system liquid tubing from the liquid container.
2. Connect the maintenance tube (30043739) to the system tubing.

3. Place the open end of the maintenance tubing in a bottle with cleaning agent.
4. Flush with cleaning agent (20 ml with RapidWash and 10 ml with diluter).
5. Leave to soak for 20 minutes.
6. Place the tubing in a bottle with DI water.
7. Rinse twice with DI water (20 ml with RapidWash and 10 ml with diluter).

WARNING

Flammable liquids!

Fire hazard caused by flammable liquids or system liquid.

- Avoid the formation and accumulation of flammable vapors.
- Do not operate the system without deck trays.

8. Place the tubing in a bottle with alcohol.
9. Flush with alcohol (20 ml with RapidWash and 10 ml with diluter).
10. Remove the maintenance tubing from the system tubing and connect the system tubing to the system liquid container.
11. Flush twice with DI water (20 ml with RapidWash and 5 times the diluter volume).
12. Check for bubbles in the tubing.
13. Flush again if bubbles are visible.

7.5.16 Connecting the System Liquid Container and Waste Container

To prepare the system liquid container and waste container, proceed as follows:




A trouble-free operation is only guaranteed if the original containers with controlling system from Tecan are used.

Prior to first use, the system liquid container must be manually rinsed out thoroughly to remove any solid debris from inside the bottle. Refer to section [“Cleaning the System Liquid Container and Waste Container”](#) [▶ 133].

- ✓ Tecan container with a capacity of more than 20 liters

1. Ensure that the liquid detection system (A, B) is connected correctly.

2. Ensure that the tubes (C, D) are connected correctly.



7.5.17 Cleaning the System Liquid Container and Waste Container

To clean the liquid container and the waste container, proceed as follows:

1. Empty the wash liquid container manually.
2. Clean the liquid container in a basin with the cleaning agent and rinse.
3. Disinfect the liquid container with alcohol.
4. Connect the system liquid and the waste container, refer to section ["Connecting the System Liquid Container and Waste Container" \[132\]](#).

7.5.18 Checking Tightness of Syringes

To check the correct tightness of the syringes proceed as follows:




Fig. 56: Checking tightness

A Plunger lock screw

B Syringe screw


1. Slightly turn the plunger lock screw (A) counterclockwise to loosen it.
2. Tighten the syringe screw (B)—i.e., turn it right.
3. Turn the plunger lock screw clockwise to tighten it.

7.5.19 Checking Gaskets (MCA)

To check gaskets (MCA), proceed as follows:

1. Inspect the gaskets (A) for damage.

Ensure that the gaskets are not damaged. A damaged gasket must be replaced. Refer to section [“Replacing Gaskets \(MCA\)”](#) [135].



7.5.20 Replacing Gaskets (MCA)





Fig. 57: Parts and tools


- | | | | |
|----------|---------------------------|----------|-----------------------|
| A | Gasket | B | Blunt tube |
| C | Blunt tube insertion tool | D | Gasket insertion tool |
| E | Blunt tube removal tool | | |

To replace gaskets (MCA), proceed as follows:

1. Remove the adapter plate.
2. Move the pipetting head to the front and raise as high as possible.
3. Push the blunt tube removal tool approx. 2 mm into the channel.
4. Secure the blunt tube by turning the tool clockwise.



5. Pull the blunt tube out of the channel.




6. Use long nose pliers to remove the gasket from the channel.




7. Lubricate the new gasket with a film of mineral oil.

8. Insert the lubricated gasket into the channel using the gasket insertion tool.



9. Lubricate the outer surface of a new blunt tube with mineral oil.
10. Insert the lubricated blunt tube gently into the channel using the blunt tube insertion tool.

The blunt tube must be fully inserted inside the channel.



11. Inspect the new gasket to ensure evenness with the other gaskets.

12. Run the **MCA leakage test** method.

7.5.21 Cleaning Fixed Tip Block (MCA)

1. Clean the fixed tip block (MCA) with alcohol, using a lint-free cloth.
2. Wipe the tip block dry with a lint-free tissue or blow the tip block dry by means of oil-free compressed air.
3. Ensure that the fixed tip block (MCA) is clean and free of deposits.
4. Store the tip block in the tip block box.

To ensure that contamination is avoided when handling the tip block:

- The tip block must be stored in a dust-free place.
- Never touch the tips with your fingers. The tip block must always be held by the PEEK block during handling.
- Never place the tip block with its tips on the table.

7.5.22 Cleaning Plate Adapter (MCA)


To clean the plate adapter, proceed as follows:

1. Check for dirt in the channels.
2. Use alcohol to remove any grease.
3. Dry with compressed air.

7.5.23 Cleaning Arm Guide

To clean the arm guide, proceed as follows:


1. Clean the arm guide roller (A) on the arm guide with a cotton swab or a lint-free cloth on a screwdriver.
2. Clean the arm rails (B) with a lint-free cloth.
3. Where present, clean the upper surface of the guide rail on the MCA arm guide with a lint-free cloth.



7.5.24 Tightening a DiTi Cone

To tighten the FCA DiTi Cone, proceed as follows:

1. Hold the tip adapter (D) and tip ejector tube (C).
2. Tighten the DiTi cone (A) using the DiTi cone wrench (B).



3. Run the **FCA Routine Maintenance** method.

7.5.25 Frida Reader

Insert




Fig. 58: Frida Reader insert

For installation fit the insert into the Frida Reader and align the markings.

Blind Plug



Fig. 59: Frida Reader blind plug

The blind plug protects the Frida Reader when the insert is removed. For installation fit the blind plug into the Frida Reader.

8 Troubleshooting

Consult this chapter for help on resuming operation after a problem has occurred with the Fluent. For further information or, in the event of problems not covered in this manual, or in insufficient detail, please consult section “Customer Support” [▶ 185].

8.1 Safety Instructions for This Chapter

CAUTION

Cross contamination due to damaged tips after crash!

Bent tips or damaged tip coating cause pipetting inaccuracy and liquid detection errors.

- Check the fixed tips after a crash. Refer to section “Checking Fixed Tips” [▶ 161].

8.2 Troubleshooting Tables



The troubleshooting tables list possible problems, causes, and corrective measures. For further information or, in the event of problems not covered in this manual, or in insufficient detail, please consult section “Customer Support” [▶ 185].

8.2.1 Instrument Troubleshooting

Tab. 28: Instrument troubleshooting table

Problem/Error	Possible Cause	Corrective Measure
System liquid leak	Tubing and/or tubing connections are not tight. Syringe is leaking.	Please consult section “Customer Support” [▶ 185].
Communication error	Power is not ON. Power or communication is interrupted. No communication.	Switch off the instrument. Wait until the instrument status lamp and power supply lamp switch off. Switch off the PC. Check cable and plugs. Switch on the instrument and the PC.
	X-, Y- or Z-drive blocked.	Check for obstacles. NOTICE! Ensure that the arms can move freely.

Problem/Error	Possible Cause	Corrective Measure
Initialization error	Arms cannot be initialized.	Check for obstacles. NOTICE! Ensure that the arms can move freely.
	Hardware problem.	Please consult section “Customer Support” [▶ 185].
Front safety panel door sensor and door lock are damaged	Mechanical failure of the door locks.	Switch off the instrument. Please consult section “Customer Support” [▶ 185].
Safety panel missing or damaged	Safety cannot be guaranteed.	Switch off the instrument. Please consult section “Customer Support” [▶ 185].
Liquid detection (cLLD) fault	Dirty contact surface. Inadequate contact between labware and segment.	Prepare deck. Refer to section Check before Starting a Method. Clean the contact surface. Refer to section “Cleaning Runners and Segments” [▶ 123].
	System Liquid has a conductivity > 10 µS/cm for cLLD compatibility.	Contact the key operator.

8.2.2 Flexible Channel Arm (FCA) Troubleshooting

Tab. 29: Flexible Channel Arm troubleshooting table

Problem/Error	Possible Cause	Corrective Measure
Loose DiTi cone CAUTION! Inaccurate pipetting volumes!	Insufficiently tightened DiTi cone.	Tighten the DiTi cone.
Disposable tip not fetched	Insufficiently tightened DiTi cone.	Tighten the DiTi cone.
Disposable tip not discarded	Insufficiently tightened DiTi cone.	Tighten the DiTi cone.
	Reused DiTis	Ensure tips are new. DiTis are not recommended for reuse.

Problem/Error	Possible Cause	Corrective Measure
Tips not aligned with labware on a single carrier	Carrier in wrong position. Segment not locked in place. Labware not positioned correctly.	Ensure correct carrier position. Refer to section “Loading Standard Runners” [▶ 78]. Lock segment in place. Refer to section “Checking Segment” [▶ 178].
Tips not aligned with labware on several carriers	Defective arm alignment caused by a collision.	Please consult section “Customer Support” [▶ 185].
Tip collides with bottom of labware	Wrong labware. Labware not positioned correctly.	Ensure that the labware on the deck corresponds to the method deck layout.
DiTi drips	Dirty DiTi cone causing leak.	Clean the DiTi cone.
	Reused DiTis	Ensure tips are new. DiTis are not recommended for reuse.
Error message: Pressure out of range (Air FCA)	Wet inline filter after aspiration with wrong DiTi size.	Ensure that the size of the DiTi on the deck corresponds to that defined in the method. Check the inline filter. Refer to section “Checking Inline Filter (Air FCA)” [▶ 154].
Error messages: DiTi not fetched DiTi not dropped	Magnetic field interfering with DiTi-presence sensor.	—

8.2.3 Mix and Pierce

Tab. 30: Troubleshooting

Symptom	Possible Cause	Corrective Measures
Piercing tip cannot be retracted with software commands	Stuck piercing tip	Refer to section “Retracting Stuck Piercing Tips” [▶ 172].
Damaged piercing tip	Bent piercing tip Damaged tip	Replace the piercing tip. Refer to sections “Removing Piercing Tips” [▶ 165] and “Installing Piercing Tips” [▶ 168].

Symptom	Possible Cause	Corrective Measures
Piercing errors	Piercing tip too dry	Lubricating with water (wash station)
	Damaged tip	Replace the piercing tip. Refer to sections “Removing Piercing Tips” [▶ 165] and “Installing Piercing Tips” [▶ 168].
	Bent piercing tip	
	Wrong piercing parameter	Contact the key operator.
	Wrong movement type used	Contact the key operator.
	Wrong tubes used	Use supported tubes. Refer to section “Tube Rotator Runners” [▶ 58].
	Arm has reached its lifetime	Please consult section “Customer Support” [▶ 185].
Liquid handling problems	Clogged piercing tips	Flush piercing tips. Check washing procedure in general.
	Damaged tip	Replace the piercing tip. Refer to sections “Removing Piercing Tips” [▶ 165] and “Installing Piercing Tips” [▶ 168].
	Syringes not properly mounted	Check tightness of syringes. Refer to section “Checking Tightness of Syringes” [▶ 133].
	Bubbles in Liquid System	Flush. Please consult section “Customer Support” [▶ 185].
Hemolysis problems	Sample dilution	Higher excess volume or partitioning volume
		Saline 0.9% as partitioning volume
		Lower pipetting speeds
	Damaged tip	Replace the piercing tip. Refer to sections “Removing Piercing Tips” [▶ 165] and “Installing Piercing Tips” [▶ 168].
	Mixing parameters	Make sure that the used rotation / oscillation parameters of the tube rotator do not lead to hemolysis

Symptom	Possible Cause	Corrective Measures
Sample in syringe	Any	Clean system. Refer to section “Cleaning Liquid Path” [▶ 131].
	Wrong air gap.	Validate wash procedure.
	Syringes not properly mounted.	Check tightness of syringes. Refer to section “Checking Tightness of Syringes” [▶ 133].
		Run the FCA Leakage Method .
		Bigger leading air gap.
		Slower aspiration speeds.
	Incorrect sample preparation for liquid handling. Sample source tubes contain solid particles like clots, cell debris, etc.	Ensure proper sample preparation to allow pipetting of sample liquid.
		Make sure that sample source tubes do not contain solid particles like clots, cell debris, etc.
	Incorrect sample preparation for liquid handling. Tubes are not properly filled and still contain partial vacuum that reduces the leading air gap during piercing.	Make sure that the sample source tubes are properly filled with the tube’s target volume.
		Make sure that sample source tubes do not contain vacuum.
		Increase the leading air gap to compensate for possible remaining vacuum.
Piercing tip bends during wash procedure	Piercing tip is not centered in wash station cleaner holes	Create a copy of the wash station and teach the pipetting positions.

Symptom	Possible Cause	Corrective Measures
Incorrect liquid level detection: only on specific channels	Bent piercing tip: Piercing Tip is bent and therefore touches the tube wall during piercing	Replace the piercing tip. Refer to sections “Removing Piercing Tips” [▶ 165] and “Installing Piercing Tips” [▶ 168].
	Piercing position is incorrect and therefore the piercing tip touches the tube wall during piercing.	Use Tecan manufactured labware. Refer to section “Tube Rotator Runners” [▶ 58]. Teach/adjust labware pipetting position
	Orientation of piercing tip is incorrect.	Mount piercing tip with opening facing to the front of the instrument. Refer to section “Installing Piercing Tips” [▶ 168].
Incorrect liquid level detection: constant deviation of expected liquid level and detected liquid level	Manufacturing tolerances of the arm, Tube Rotator and/or instrument in combination with piercing force can lead to noticeable Z-offset during liquid level detection.	Teach/adjust the custom attribute "PiercingDetection-HeightCompensation" in the tube labware definition

8.2.4 Multiple Channel Arm (MCA) Troubleshooting

Tab. 31: Multiple Channel Arm troubleshooting table

Problem/Error	Possible Cause	Corrective Measure
Tips are not aligned with the carriers	Mechanical fault	Please consult section “Customer Support” [▶ 185].
	Arm crash	
Microplate and pipetting head not 100% parallel	Crash	Contact the key operator to check parallelism of pipetting head and the deck.
During pipetting, the pipetting head stops generating an error	Aspiration and dispensing acceleration is too fast compared to speed.	Acceleration must be in a reasonable relation to aspirating and dispensing speed.
	Aspiration and dispensing deceleration is too fast compared to speed.	Deceleration must be in a reasonable relation to aspirating and dispensing speed. Problem cannot be solved. Please consult section “Customer Support” [▶ 185].

Problem/Error	Possible Cause	Corrective Measure
Several or all pipetting channels leaking	Wrong disposable tips, fixed tip block, fixed tip adapter, tip cone seals or gaskets.	Always use disposable tips, fixed tip block, fixed tip adapter, tip cone seals or gaskets supplied by Tecan. Contact the key operator to flush the system and check the system for leaks.
	Tip cone seals or gaskets old or defective.	Please consult section “Customer Support” [▶ 185]. Contact the key operator to flush the system and check the system for leaks.
	Pipetting head is faulty.	Please consult section “Customer Support” [▶ 185].
Single channel leaking	Tip cone seal, gaskets or other seals in the pipetting head are defective.	Contact the key operator to flush the system and check the system for leaks. Please consult section “Customer Support” [▶ 185].
Single disposable tip not picked up correctly	Individual disposable tip is defective. Tip cone seal on this disposable tip position is defective.	Replace disposable tips. Problem cannot be solved. Please consult section “Customer Support” [▶ 185].
Single disposable tip not dropped	Individual disposable tip is defective. Tip cone seal on this disposable tip position is defective.	Problem cannot be solved. Please consult section “Customer Support” [▶ 185].
	Wrong humidity	Ensure the humidity is within operating humidity limits. Refer to section “Environmental Conditions” [▶ 44].
Several or all disposable tips not dropped	Wrong disposable tips used.	Always use disposable tips supplied by Tecan. Problem cannot be solved. Please consult section “Customer Support” [▶ 185].

Problem/Error	Possible Cause	Corrective Measure
Disposable tip box is lifted up with disposable tips when picking up disposable tips	Carrier is not adjusted correctly. X- and/or Y-offset specified incorrectly.	Adjust all carriers (mechanical) precisely. Replace the disposable tip carrier. Problem cannot be solved. Please consult section “Customer Support” ► 185].
	Disposable tip box does not meet specifications.	Always use disposable tip boxes that comply with the Society of Biomolecular Screening standards. Problem cannot be solved. Please consult section “Customer Support” ► 185].
	Disposable tip carrier is defective (malfunction on disposable tip box retainers).	Always use disposable tip boxes that comply with the Society of Biomolecular Screening standards. Problem cannot be solved. Please consult section “Customer Support” ► 185].
Inaccurate pipetting results	Disposable tips are not picked up properly. Liquid handling parameters are incorrect. Carriers not correctly adjusted. Pipetting head is faulty.	Contact the key operator to check the application script and the carriers. Contact the key operator to check the environment parameters and the dispensing height. Check the tip coating. If necessary, replace tip block.
	Coating of standard tips is damaged.	Contact the key operator to check the application script and the carriers. Check the tip coating. If necessary, replace tip block.
Carry-over	Coating of standard tips is damaged.	Check coating. Replace fixed tip block or fixed tip adapter.
	Unsuitable application script.	Contact the key operator to adjust the wash or liquid handling settings and use a different wash buffer.
	Wash channels are clogged.	Clean wash station. Problem cannot be solved. Please consult section “Customer Support” ► 185].

8.2.5 Robotic Gripper Arm (RGA) Troubleshooting

Tab. 32: Robotic Gripper Arm troubleshooting table

Problem/Error	Possible Cause	Corrective Measure
Microplate not picked up	No microplates on carrier. Gripper fingers cannot pick up the microplate.	Place microplate on carrier. Set gripper position. Clean RGA gripper fingers.

Problem/Error	Possible Cause	Corrective Measure
Unusual noise during arm movement	Parts are damaged or worn.	Please consult section “Customer Support” [▶ 185].
Eccentric Gripper Fingers misaligned	Spare fingers crash. Finger screws not tight enough.	Align the Eccentric Gripper Fingers. Refer to section “Checking Gripper Finger Alignment” [▶ 173]. Use a torque screw driver to tighten the screws to 3 Nm, as described in section “Basic Gripper Fingers Alignment for FES Gripper Fingers” [▶ 174].

8.2.5.1 Robotic Gripper Arm with long Z-axis (RGA-Z) Troubleshooting

Tab. 33: Robotic Gripper Arm with long Z-axis troubleshooting table

Problem/Error	Possible Cause	Corrective Measure
Microplate not picked up	No microplates on carrier. Gripper fingers cannot pick up the microplate.	Place microplate on carrier. Set gripper position. Clean RGA gripper fingers.
	Gripper fingers are slippery.	Clean RGA gripper fingers.
Unusual noise during arm movement	Parts are damaged or worn.	Please consult section “Customer Support” [▶ 185].

8.2.6 Wash System Troubleshooting

Tab. 34: Wash system troubleshooting table

Problem/Error	Possible Cause	Corrective Measure
Erroneous overflow and empty errors during process run	Overflow sensor is not clean. Sensors are not connected or faulty.	Check that sensor is connected. Refer to sensor manufacturer's manual.
Regular wash station overflow	Waste pump is defective. Waste pump tubings are kinked, clogged, or damaged.	Check waste tubing. Replace waste tubing if necessary. Please consult section “Customer Support” [▶ 185].

Problem/Error	Possible Cause	Corrective Measure
No wash liquid pumped through the wash station or the wash system is empty	Waste pump tubings are kinked, clogged or damaged.	Check waste tubing. Replace the waste pump if necessary. Please consult section “Customer Support” [▶ 185].
	Wash station is not connected. Wash container(s) empty or missing.	Refill or replace the wash container(s). Connect wash station (MCA) properly. Refer to section “Connecting Wash Station (MCA)” [▶ 152].
	Waste pump is defective.	Check waste pump. Replace the waste pump if necessary. Please consult section “Customer Support” [▶ 185].
Wash station overflow	Waste tube is below liquid surface in waste container.	Use a wash container with fixed wash tubing inlet.
	Disposable tips or algae block the wash station.	Clean the wash station. Refer to section “Cleaning the System Liquid Container and Waste Container” [▶ 133].
	Waste tubing is kinked.	Check tubing for kinks. Refer to section “Checking the Tubing on System Liquid Container and Waste Container” [▶ 84].

8.2.7 Fluent ID Troubleshooting

Tab. 35: Fluent ID troubleshooting table

Problem/Error	Possible Cause	Corrective Measure
Barcode not read	Barcode label not facing the scanner.	Unload the tube runner, turn the tubes so that the barcode labels face left. Reload the tube runner on the Fluent.
	Runner loaded too fast.	Unload the tube runner and load it again slowly.
	Poor label quality.	Enter the barcode manually or report the problem to the key operator.
	Scanner window is dirty.	Clean the scanner window. Refer to section “Weekly System Care” [▶ 111].
	Reflector is dirty.	Clean the reflector. Refer to section “Weekly System Care” [▶ 111].
	Barcode type or barcode length not pre-defined for the method.	Report the problem to the key operator.

Problem/Error	Possible Cause	Corrective Measure
Tube presence not detected	Barcode label position too low on the tube.	Report the problem to the key operator.

8.2.8 Software Troubleshooting

Tab. 36: Software troubleshooting table

Problem/Error	Possible Cause	Corrective Measure
User login screen not displayed when envisaged.	User management has not been activated in FluentControl.	Contact the key operator to activate user management.
User cannot log in.	Password is incorrect or account is locked.	Contact the key operator to reset the password or account.
Not all service actions completed. Warning appears at each FluentControl startup.	Not all envisaged service actions are marked as complete in the instrument configuration.	Please consult section " Customer Support " [▶ 185].
Touchscreen does not react to touch.	Software driver not installed.	Contact computer administrator for installation of the drivers on the installation CD and configuration of the touchscreen.
	Touchscreen interface incorrectly configured.	Open the touchscreen driver settings and ensure that the touchscreen is correctly mapped.
Touch interface is not displayed on touchscreen.	Touchscreen was not on when the software was started.	Turn on the instrument and restart the software or check the Touch Tool settings in the FluentControl configuration system.
Error on FluentControl startup.	FluentControl (SystemSW.exe) is already running in the background (Task Manager).	Open the Task Manager, process <code>SystemSW.exe</code> and restart FluentControl. Or restart the computer.
FluentControl does not communicate with connected hardware devices.	FluentControl is not properly configured for communication with hardware devices.	Contact the person responsible for configuring the system to activate the I/O state of the hardware devices.

8.3 Troubleshooting Activities


8.3.1 Connecting Wash Station (MCA)

To connect the MCA Wash System wash station, proceed as follows:

1. Switch off the instrument.
2. Connect the wash tube (A) and the waste tube (B) to the wash station.



3. Connect the wash tube (A) and the waste tube (B).



4. Connect all tubes (C) to the MCA wash control unit.




5. Contact the key operator to check operation readiness.


8.3.2 Position Labware

To ensure that labware is positioned correctly in the nest for precise arm access, proceed as follows:

1. Place the labware (A) on the nest (B).



2. Push the labware gently against the static positioner (C).



3. Slide the diagonal sliding positioner (E) towards or away from the labware to precisely fit the labware.
4. Slide the vertical and horizontal sliding positioners (D, F) towards or away from the labware to fix the labware.
5. Lift the labware from the nest.

Ensure that there is no friction when placing or removing the plate.

8.3.3 Checking Inline Filter (Air FCA)

A control system is installed on each channel to protect Air FCA pipetting channels against over-aspiration of liquid.

- ✓ The **Air FCA Routine Maintenance** method includes an inline filter check, which detects wet, damaged or mispositioned filters as well as missing filters.
1. Run the **Air FCA Routine Maintenance** method to check the inline filter inside the DiTi cone of an Air FCA pipetting channel.

In the event of an error, the inline filter must be changed. Refer to section [“Changing Inline Filter \(Air FCA\)” \[154\]](#).

8.3.4 Changing Inline Filter (Air FCA)

To change the inline filter, proceed as follows:

- ✓ Decontaminated disposable tip cone.
- ✓ DiTi cone removed from the channel. Refer to section [“Removing DiTi Cone \(Air FCA\)” \[156\]](#).





Fig. 60: Removing inline filter

- A** DiTi cone
C Filter removal tool
B Inline filter

1. Pierce the inline filter (B) sideways with the filter removal tool (C).
2. Pry out the inline filter with the filter removal tool. Note that the filter may be contaminated with process liquids.
3. Dispose of the inline filter.




4. Clean the DiTi cone (A) with alcohol.
DiTi cone must be dry before refitting.
5. Place the new inline filter on a clean and flat surface.
6. Press the inline filter into the DiTi cone.
The inline filter must not protrude from the DiTi cone.
7. Check the inline filter according to the method defined by your key operator.

8.3.5 Removing DiTi Cone (Air FCA)


To remove the DiTi Cone (Air FCA), proceed as follows:

- ✓ DiTi cone wrench available.
- 1. Switch off the instrument.
- 2. Open front safety panel.
- 3. Manually raise all the Z-rods to their uppermost position.
- 4. Move all the Z-rods towards the front of the instrument.
- 5. Spread the Z-rods as wide as possible.
- 6. Hold the tip adapter (D) and tip ejector tube (C).
- 7. Unscrew the DiTi cone (A) using the DiTi cone wrench (B).



- 8. Pull down the DiTi cone carefully.


In some cases the tip ejector tube (C) or the adapter cylinder (B) may still be attached to the DiTi cone (A). Refer to section [“Assembling DiTi Ejector Tube \(Air FCA\)”](#) [157].



8.3.6 Assembling DiTi Ejector Tube (Air FCA)

To assemble the DiTi ejector tube (Air FCA), proceed as follows:

- ✓ The DiTi ejector tube has been removed, according to the instructions.
 - ✓ DiTi cone wrench available.
1. Insert the sealing sleeve (B) in the adapter cylinder (A).



2. Screw the DiTi cone (C) to the assembled cylinder.



3. Screw the tip ejector tube (D) to the assembled cylinder.




8.3.7 Installing DiTi Cone (Air FCA)

To install the Air FCA DiTi Cone, proceed as follows:

- ✓ DiTi cone is fully assembled: Refer to section [“Assembling DiTi Ejector Tube \(Air FCA\)”](#) [▶ 157].
 - ✓ DiTi cone wrench available.
1. Slide the adapter cylinder into the tip ejector tube (C).
 2. Hold the tip adapter (D) and tip ejector tube (C).

3. Screw in the DiTi cone (A) using the DiTi cone wrench (B).




4. Run the **Air FCA Routine Maintenance** method.

8.3.8 Removing the DiTi Option (FCA)

To remove the DiTi option, proceed as follows:

✓ DiTi cone wrench

1. Switch off the instrument.
2. Open the front safety panel.
3. Manually raise all the Z-rods to their uppermost position.
4. Move all the Z-rods towards the front of the instrument.
5. Spread the Z-rods as wide as possible.
6. Hold the tip adapter (D) and tip ejector tube (C).
7. Unscrew the DiTi cone (A) with the DiTi cone wrench (B).




8. Pull down the DiTi cone carefully.

8.3.9 Installing the DiTi Option (FCA)

To install the DiTi option, proceed as follows:


1. Screw the adapter cylinder (B) in the tip ejector tube (C).
2. Screw the DiTi cone (A) to the adapter cylinder. Use the DiTi cone wrench.



3. Push the tubing onto the plastic needle until the tubing attaches firmly to the DiTi option.



4. Hold the tip adapter (D) and tip ejector tube (C).
5. Screw in the DiTi cone (A) using the DiTi cone wrench (B).



8.3.10 Checking Fixed Tips

To check the fixed tips, proceed as follows:

NOTICE

Pipetting inaccuracy and liquid detection errors!

Bent or damaged tip coating causes pipetting inaccuracy and liquid detection errors.

- Never work with damaged or bent tips.

1. Switch off the instrument.
2. Open front safety panel.
3. Inspect the fixed tips.
4. Inspect the fixed tip coating with a mirror.

Ensure that the fixed tips are not bent. If the fixed tip coating is damaged or the fixed tip is bent, it must be replaced. Refer to section “Removing Fixed Tips” [► 162].

8.3.11 Removing Fixed Tips

To remove fixed tips, proceed as follows:

- ✓ Fixed tips have been cleaned. Refer to section “System Care Tables” [► 106].
- ✓ Fixed tips have been checked. Refer to section “Checking Fixed Tips” [► 161].




Fig. 61: Standard tip removal

- | | | | |
|----------|------------------|----------|-------------|
| A | Z-rod | B | Tip adapter |
| C | Lock nut | D | Tip |
| E | Pipetting tubing | | |

1. Switch off the instrument.
2. Open front safety panel.
3. Manually raise all the Z-rods (A) to their uppermost position.
4. Spread the Z-rods as wide as possible.
5. If an adjustable fixed tip is installed, loosen the four tip adjustment screws.
6. Unscrew lock nut (C) while holding the fixed tip immediately below the lock nut with the other hand.
7. Remove the lock nut (C), sliding it along the tip axis.
Avoid contact between the lock nut and the tip coating.
8. If the tip (D) is adjustable, turn the lock nut (C) upside down on a clean surface, and remove the O-ring and the washer.

9. If the channel is equipped with the low volume option, unscrew the flange on top of the solenoid valve to free the pipetting tubing (E) running through the Z-rod (A).
10. Extract the pipetting tubing (E) a certain distance (a) out of the tip adapter (B) by pulling on the tip (D).

Use a dry emery cloth for improved grip on the pipetting tubing—not on the tip.

8.3.12 Installing Fixed Tips




Fig. 62: Standard tip installation

- | | | | |
|----------|------------------|----------|-------------|
| A | Z-rod | B | Tip adapter |
| C | Pipetting tubing | D | Tip |
| E | Lock nut | | |

To install fixed tips, proceed as follows:

1. Carefully pull the pipetting tubing approx. 25 mm (1 in.) (a) out of the tip adapter.

Use a small piece of emery cloth to grip the tubing near the end to ensure a better grip.

If a tip has been installed before, cut approx. 5 mm (0.2 in.) (b) from the pipetting tubing, using a sharp knife to obtain a straight cut.




Fig. 63: Te-PS tubing widener

C	Pipetting tubing	D	Te-PS tip
F	Emery cloth	G	Te-PS tubing widener
H	Washer, white (FEP)	I	O-ring, black
J	Adjustable lock nut	K	Tip adjustment screw

2. In case of Te-PS tips or low-volume tips:

Use the Te-PS tubing widener (G) to widen the tubing end by pushing the Te-PS tubing widener up to the hilt into the tubing while turning the tool.

While the tubing is still wide, push the Te-PS tip into the tubing end by approx. 4 mm (0.16 in.).

3. Slide the lock nut onto the tip.

If the tip is adjustable (e.g., Te-PS), slide the lock nut over the washer (H) and the O-ring (I).

NOTICE! Avoid contact with the delicate end of the tip and its coating.

4. Insert the tip and the pipetting tubing in the tip adapter.

5. Screw the lock nut onto the tip adapter and tighten.

If the tip is adjustable (e.g., Te-PS), tighten the lock nut so that the four tip adjustment screws (K) are at a 45° angle to the deck's X/Y coordinate system.

6. Clean the fixed tips. Refer to section [“End of Day”](#) [108].

7. Run a pipetting precision test as defined by the key operator.

8.3.13 Removing Piercing Tips




Fig. 64: Piercing tip removal

A	Z-rod	B	Tip adapter
C	Lock nut	D	Piercing tip
E	Piercing tip protection	F	Pipetting tubing
x	25 mm (1 in.)		

To remove piercing tip, proceed as follows:

- ✓ The instrument is switched off.
- 1. Open the front safety panel.
- 2. Manually raise all the Z-rods to their uppermost position.
- 3. Move all the Z-rods towards the front of the instrument.
- 4. Spread the Z-rods as wide as possible.

5. Cover the piercing tips with piercing tip protections. Start with the rearmost piercing tip.




6. Unscrew the lock nut, holding the piercing tip immediately below the lock nut with the other hand.

7. Pull the pipetting tubing approx. 25 mm out of the tip adapter by pulling on the tip. Hold the piercing tip at its upper end when pulling.



8. Pull the piercing tip off the tubing, withholding the tubing with the other hand.
9. Do not remove the piercing tip protection. Discard it with the piercing tip into the biological waste container.



8.3.14 Installing Piercing Tips





Fig. 65: Piercing tip installation

- | | | | |
|----------|------------------|----------|----------------|
| A | Z-rod | B | Tip adapter |
| C | Pipetting tubing | D | Lock nut |
| E | Piercing tip | F | Tip protection |
| x | 25 mm (1 in.) | | |

To install piercing tips, proceed as follows:


- ✓ The instrument is switched off.
 - ✓ A key operator is available.
1. Open the front safety panel.
 2. Manually raise all the Z-rods to their uppermost position.
 3. Move all the Z-rods towards the front of the instrument.
 4. Spread the Z-rods as wide as possible.

5. Open the piercing tip packaging.
Do not remove the tip protection (F).




Installation order of the piercing tips: from back to front

6. Carefully pull the pipetting tubing approx. 25 mm out of the tip adapter.
7. Push the blank, conical end of the piercing tip into the tubing end.




8. Insert the piercing tip and the pipetting tubing in the tip adapter.


9. Screw the lock nut onto the tip adapter and tighten it by hand.



10. Open the lock nut slightly. Move the tip protection slightly down in order to access the piercing tip shaft. Do not completely remove the tip protection yet.
11. Turn the piercing tip until the tip opening faces the front side of the instrument. Keep the piercing tip in this orientation with one hand and tighten the lock nut with the other hand.
12. Check that all tip openings face the front side of the instrument.



13. After installing all piercing tips remove all tip protection. Start with the rearmost piercing tip.



14. Contact a key operator for resetting the counter in FluentControl.
15. Contact a key operator for performing a QC kit test. Refer to [“Reference Documents” \[12\]](#).
16. Run the **Piercing FCA Leakage** method.
17. Run a pipetting precision test (recommendation: use QC kit) as defined by the key operator.


8.3.15 Retracting Stuck Piercing Tips




In case a piercing tip gets stuck so that it cannot be retracted with software commands, it has to be removed manually.

To retract stuck piercing tips, proceed as follows:

- ✓ The instrument is switched off.
- 1. Open the front safety panel.
- 2. Manually raise all retracted Z-rod to their uppermost position.
- 3. Cover all retracted piercing tips with piercing tip protections (A). Start with the rearmost piercing tip.



4. Place the piercing tip removal tool next to the stuck tip on a robust and stable surface and fit it under the lock nut.



5. Turn the knob of the piercing tip removal tool until the tip is completely retracted.
6. Turn the knob in the opposite direction and lower the retractor about 1 cm.
7. Remove the piercing tip removal tool.
The stuck piercing tip is now retracted.
8. Clean the piercing tip removal tool with alcohol.
9. Check the piercing tip for any damage (e.g., bent piercing tip, damaged tip).
10. Replace the piercing tip if it is damaged. Refer to section “[Removing Piercing Tips](#)” [▶ 165] and section “[Installing Piercing Tips](#)” [▶ 168].
11. Remove all piercing tip protections by holding the lock nut with one hand and removing the tip protections with the other hand. Start with the rearmost piercing tip.
12. Clean the piercing tip removal tool with alcohol.

Also see about this

- ▢ “[Removing Piercing Tips](#)” ▶ 165
- ▢ “[Installing Piercing Tips](#)” ▶ 168

8.3.16 Checking Gripper Finger Alignment

Realignment of gripper fingers may be needed after a crash or when implementing spare Gripper Fingers.



Misalignment after a crash:


- Analyze the situation.
- Evaluate potential causes of the crash or cause of the finger misalignment such as misaligned drawer of a reader, washer, incorrectly taught/positioned hotel or another segment.
- Select a procedure below on the basis of precision requirements.

1. If Gripper Fingers do not need to satisfy above-average precision, perform a basic alignment. Refer to section “[Basic Gripper Fingers Alignment for FES Gripper Fingers](#)” [▶ 174], or to section “[Basic Gripper Fingers Alignment for Fixed Gripper Fingers](#)” [▶ 175].

2. If Gripper Fingers must satisfy advanced requirements (Z-deviation $< \pm 0.2$ mm), perform the Advanced Gripper Finger alignment procedure. Refer to section [“Advanced Gripper Fingers Alignment for FES Gripper Fingers”](#) [▶ 176] or [“Advanced Gripper Fingers Alignment for Fixed Gripper Fingers”](#) [▶ 176].
3. The Gripper fingers can be mounted with two different screws:
 - a) Torx screw M4x12, tightened with a torque of 3 Nm.



- b) Allen screw M4x12 in combination with a tension washer (observe position according to the illustration below), tightened with a torque of 3.5 Nm.




If no torque screwdriver is available, tighten screw until the washer is pressed flat and resistance increases. Then tighten $\frac{1}{12}$ rotation further. This corresponds to approx. 3.5 Nm.

8.3.17 Basic Gripper Fingers Alignment for FES Gripper Fingers

For a basic alignment, proceed as follows:

- ✓ Misalignment is clearly visible.
 - ✓ No above-average precision is required.
 - ✓ Torque screwdriver available (with a 3 or a 3.5 Nm option).
If no torque screwdriver is available:
Torx screw: tighten the screws firmly but do not use excess force.
Allen screw: refer to [“Checking Gripper Finger Alignment”](#) [▶ 173].
1. Remove Gripper Finger from the gripper head.
 2. Loosen the screw between Gripper Finger and FES finger adapter.

3. Press the Gripper Finger against the upper and the rear stop of the adapter as shown in the illustration below and tighten the screw with a torque screwdriver (3 or 3.5 Nm).




8.3.18 Basic Gripper Fingers Alignment for Fixed Gripper Fingers

For a basic alignment, proceed as follows:

- ✓ Misalignment is clearly visible.
- ✓ No above-average precision is required.
- ✓ Torque screwdriver available (with a 3 or 3.5 Nm option).


1. Loosen the fixing screw (A).



2. Press the Gripper Finger against the upper and the rear stop of the Gripper Head Mount and tighten the screw with a torque screw driver (3 or 3.5 Nm).

8.3.19 Advanced Gripper Fingers Alignment for FES Gripper Fingers

- ✓ Torque screwdriver available (with a 3 or 3.5 Nm option).
If no torque screwdriver is available:
Torx screw: tighten the screws firmly but do not use excess force.
Allen screw: refer to [“Checking Gripper Finger Alignment”](#) [▶ 173].
- 1. Use the Move Tool to move the Z-height to a height of about 3 mm above the worktable.
- 2. Loosen the screw between the Gripper Finger and the FES finger adapter.



3. Ensure that the FES finger adapter is firmly connected to the gripper head.
The finger adapters are held in place in one direction by a magnet.
4. Use the Move Tool to move the Z-height to a height of 0 mm above the worktable.
Move slowly for the last tenth of a millimeter.
Note: If you do not know how to access or operate the Move Tool, contact your key operator.
5. Press the Gripper Finger against the gripper head and the reference surface as shown in the illustration and tighten the screw with a torque of 3 or 3.5 Nm.
6. Check the adjustment by rotating the head to 90°, 180°, 270°, by hand. A misalignment at the different positions indicates a misalignment of the head or the arm. In this case, an FSE needs to check the alignment.

8.3.20 Advanced Gripper Fingers Alignment for Fixed Gripper Fingers

- ✓ Torque screwdriver available (with a 3 or 3.5 Nm option).
- 1. Use the Move Tool to move the Z-height to a height of about 3 mm.
- 2. Loosen the screw between the Gripper Finger and the gripper head.
- 3. Use the Move Tool to move the Z-height to a height of 0 mm.
Move slowly for the last tenth of a millimeter.
Note: If you do not know how to access or operate the Move Tool, contact your key operator.

4. Press the Gripper Finger against the gripper head and the reference surface and tighten screw with a torque of 3 or 3.5 Nm.
5. Check the adjustment by rotating the head to 90°, 180°, 270°, by hand. A misalignment at the different positions indicates a misalignment of the head or the arm. In this case an FSE needs to check the alignment.

8.3.21 Checking Segment

Check that segment is closed.




Fig. 66: Segment closed




Fig. 67: Segment open

8.3.22 Removing Positioning Pins

To remove positioning pins, proceed as follows:

NOTICE

Crash or process error!

Crash and other process error can be a result of inaccurate positioning of elements on a deck segment due to loose positioning pins.


- Do not operate the Fluent when positioning pins are missing.

✓ Positioning pins are broken.

1. Slide the pin remover over positioning pin.



2. Lift the pin remover handle and pull the pin out of the deck segment.



8.3.23 Replacing Lock Pins and Positioning Pins

NOTICE

Crash or process error!


Crash and other process error can be a result of inaccurate positioning of elements on a deck segment due to loose positioning pins.

- Do not operate the Fluent when positioning pins are missing.

✓ Positioning pins indicated in the system care table are available.

1. Remove segment from the deck:
Refer to section [“Removing Segments” \[76\]](#).
2. Press the new lock pin (B) in the hole (A).

3. Press the new positioning pin (C) in the hole (A).



4. Place segment on the deck:
Refer to section [“Placing Segments”](#) [► 75].

9 Packing, Unpacking, Transport, Storage and Disposal

This chapter includes regulatory information about recycling and packaging labels that must be followed.

NOTICE

Prevent damage by unqualified and unauthorized personnel!

Packing, unpacking, transport and storage may only be performed by Tecan personnel or personnel authorized by Tecan!

- Please consult the [“Customer Support”](#) [► 185].

For information about moving the instrument, refer to section [“Moving the Instrument on a Cabinet within the Laboratory”](#) [► 114].

9.1 Packaging Labels

Correct and complete marking of packaging helps to prevent incorrect handling, accidents, incorrect delivery, loss of weight and damage during storage.

Tab. 37: Packaging symbols

Symbol	Meaning	Description
	Recycle	The packaging material can be recycled. Do not dispose of as domestic waste. Information on the material used for this packaging is provided beneath the symbol.
	This side up	Ensure that the package is transported and stored with the top side, indicated by the arrows, uppermost. Do not topple over.
	Keep dry	Ensure that the package does not get wet during transport and storage.
	Fragile	Handle the package with care. There are fragile goods inside.
	Keep away from sunlight	Ensure that the package will not be exposed to heat during transport and storage. Protect against strong sunlight.
	Do not stack	Do not stack packages. The package is not designed to carry extra weight.

9.2 Disposal

This section includes regulatory information about recycling that must be followed.

NOTICE


Recycling in accordance with applicable legal regulations!

Observe the laws applicable in your country for recycling.

9.2.1 Local Requirements European Union

The European Commission has released the Directive on Waste Electrical and Electronic Equipment (WEEE; 2012/19/EU).

Since August 2005, producers have been responsible for taking back and recycling electrical and electronic equipment.

Marking	Explanation
	Negative environmental impacts associated with the treatment of waste. - Do not treat electrical and electronic equipment as unsorted municipal waste. - Collect waste electrical and electronic equipment separately.

9.2.2 Local Requirements People's Republic of China

Marking for the Restriction of the Use of Hazardous Substances in Electronic and Electrical Products

The People's Republic of China Electronic Industry Standard SJ/T11364-2014 **Marking for the Restriction of the Use of Hazardous Substances in Electronic and Electrical Products** requires the marking for the restriction of the use of hazardous substances in electronic and electrical products.

In accordance with the requirements specified in SJ/T11364-2014, all electronic and electrical Tecan products sold in the People's Republic of China are labeled with a marking for the restriction of the use of hazardous substances.

Marking	Explanation
	This marking indicates that this electronic product contains certain hazardous substances and can be safely used during the environment-friendly use period, but it shall enter the recycling system after the environment-friendly use period.

9.2.3 Other Requirements

Marking	Explanation
	This lamp contains mercury • Recycle or dispose of as required by applicable local laws.

10 Customer Support

This paragraph explains which files and information Tecan requires to perform a first assessment of an issue.

If you have any comments on this Operating Manual or suggestions for improvement, please send them by e-mail to docfeedback@tecan.com. In your e-mail, please specify the manual name, the document ID and the manual version. This information is shown at the bottom of each printed page and on the first page of the help file (context-sensitive help of software products).

10.1 Contacts

Please contact your local distributor or one of the addresses below.

Also see our homepage on the web: www.tecan.com

Tab. 38: Customer Support contacts

Country/Region	Address	Telephone/Telefax/E-mail	
Asia	Tecan Asia Pte Ltd. 18 Boon Lay Way, #10-106 TradeHub 21 Singapore 609966 Singapore	Phone Fax E-mail	+65 6444 1886 +65 6444 1836 tecan@tecan.com.sg
Australia New Zealand Pacific Islands	Tecan Australia Pty Ltd 21 / 3 Westside Avenue Port Melbourne Vic 3207 Australia	Phone Phone Fax E-mail	Toll Free: 1300 808 403 +61 3 9647 4100 +61 3 9647 4199 helpdesk-aus@tecan.com
Austria	Tecan Austria GmbH Untersbergstrasse 1a 5082 Grödig Austria	Phone Fax E-mail	+43 6246 8933 256 +43 6246 72770 helpdesk-at@tecan.com
Belgium	Tecan Benelux B.V.B.A. Mechelen Campus Schaliënhoevedreef 20A 2800 Mechelen Belgium	Phone Fax E-mail	+32 15 42 13 19 +32 15 42 16 12 tecan-be@tecan.com
China	Tecan (Shanghai) Trading Co., Ltd. Room 1802, 1803, 1804 and Room 205, HongJia Tower, 388 Fushan Road, Pudong New Area, Shanghai, P.R.China	Phone Fax E-mail	+86 21 2206 32 06 +86 40 0821 38 88 +86 21 2206 52 60 helpdesk-cn@tecan.com

Country/Region	Address	Telephone/Telefax/E-mail	
France	Tecan France S.A.S.U Tour Swiss Life 1 bd Marius Vivier Merle F- 69 003 Lyon France	Phone Fax E-mail	+33 4 72 76 04 80 +33 4 72 76 04 99 helpdesk-fr@tecan.com
Germany	Tecan Deutschland GmbH Werner-von-Siemens-Straße 23 74564 Crailsheim Germany	Phone Fax E-mail	+49 1805 8322 633 or +49 1805 TECAN DE +49 7951 9417 92 helpdesk-de@tecan.com
Italy	Tecan Italia, S.r.l. Via Brescia, 39 20063 Cernusco Sul Naviglio (MI) Italy	Phone Fax E-mail	+39 800 11 22 91 +39 (02) 92 72 90 47 helpdesk-it@tecan.com
Japan	Tecan Japan Co., Ltd. Kawasaki Tech Center 580-16, Horikawa-cho, Saiwai-ku Kawasaki, Kanagawa 212-0013 Japan	Phone Fax Phone E-mail	+81 44 556 7311 (Kawasaki) +81 44 556 7312 (Kawasaki) +81(0) 6305 8511 (Osaka) helpdesk-jp@tecan.com
Netherlands	Tecan Benelux B.V.B.A. Industrieweg 30 NL-4283 GZ Giessen Netherlands	Phone Fax E-mail	+31 20 708 4773 +31 183 44 80 67 helpdesk.benelux@tecan.com
Scandinavia	Tecan Nordic AB Sveavägen 159, 1tr SE-113 46 Stockholm Sweden	Phone Fax E-mail	+46 8 750 39 40 +46 8 750 39 56 info@tecan.se
Spain Portugal	Tecan Ibérica Instrumentación S.L. C/ Lepanto 151 Bajos E-08013 Barcelona Spain	Phone E-mail	+34 93 595 25 31 helpdesk-sp@tecan.com
Switzerland	Tecan Schweiz AG Seestrasse 103 8708 Männedorf Switzerland	Phone Fax E-mail	+41 44 922 82 82 +41 44 922 89 23 helpdesk-ch@tecan.com

Country/Region	Address	Telephone/Telefax/E-mail	
United Kingdom	Tecan UK Ltd. Theale Court 11-13 High Street Theale, Reading, RG7 5AH United Kingdom	Phone Fax E-mail	+44 118 930 0300 +44 118 930 5671 helpdesk-uk@tecan.com
USA	Tecan US, Inc. 9401 Globe Center Drive, Suite 140, Morrisville, NC 27560 USA	Phone Fax Phone E-mail	+1 919 361 5200 +1 919 361 5201 Toll Free in the US: +1 800 TECAN US or +1 800 832 2687 helpdesk-us@tecan.com
USA (Tecan Systems)	Tecan Systems, Inc. 2450 Zanker Road San Jose, CA 95131 USA	Phone Fax E-mail	+1 408 953 3100 Toll Free: +1 800 231 0711 +1 408 953 3101 helpdesk-sy@tecan.com

Abbreviations

ADT

Air Displacement Technology

Air FCA

Flexible Channel Arm with air system

ASM

Application Software Manual

CE

Conformité Européenne

cLLD

Capacitive Liquid Level Detection

CNS

Common Notification System

CSA

Canadian Standard Association

DiTi

Disposable Tip

EMC

Electromagnetic Compatibility

EN

European Norm

FCA

Flexible Channel Arm

FES

Finger Exchange System

FSE

Field Service Engineer

GLP

Good Laboratory Practice

HEPA

High-Efficiency Particulate Ar-
restance

IEC

International Electrotechnical Com-
mission

IQ

Installation Qualification

ISO

International Organization for Stan-
dardization

LED

Light Emitting Diode

Liquid FCA

Flexible Channel Arm with liquid sys-
tem

MCA

Multiple Channel Arm

MCH

Multiple Channel Head

MIO

Monitored Incubators Option

MP

Microplate

OM

Operating Manual

OQ

Operating Qualification

PC

Personal Computer

PP

Polypropylene

rcf

relative centrifugal force

RF

Radio Frequency

RGA

Robotic Gripper Arm

RGA long Z

Robotic Gripper Arm long height

RGA standard Z

Robotic Gripper Arm standard height

RUO

Research Use Only

RWP

RapidWash Pump

SN

Serial Number

Te-Shake

Tecan Shaker

Te-VacS

Tecan Vacuum Separator

USB

Universal Serial Bus

WEEE

Waste Electrical and Electronic
Equipment

WRC

Wash and Refill Center